ANYRIDGE® by MEGA'GEN





What is the AnyRidge way?

For clinicians...

less invasive, fast, simple, predictable, & esthetically superior implant treatment

Realising the ONE-DAY Implant™

For patients...

strong new esthetic & functional teeth via painless & rapid treatment

AnyRidge does it right!

AnyRidge goes FAR BEYOND standard expectations of dental implants...

The key benefits of AnyRidge implants become evident when considering immediate loading...

With the new loading protocol developed based on clinical results with AnyRidge, your patients have new smiles faster than ever....

- Guaranteed excellent stability, even with compromised bone density
- Less reduction & more preservation of cortical bone
- Wider implant possibilities than crestal width
- Clinically proven safety
- Faster & stronger osseointegration
- Esthetic design & varied abutment selection
- Super implant-prosthetic connection
- Innovative R2GATE software for completing implant
 & prosthesis in ONE DAY



AnyRidge – a new design standard on the global stage

Launched in 2009 as a biologically-inspired implant concept, AnyRidge consistently surpasses clinical benchmarks





Key Advantages

Excellent initial stability even at the compromised bone density

No screw loosening guaranteed!

Unique and valuable ISQ pattern; essential for predictable immediate or early loading.

- 004 Characteristics & Advantages
- 004 **I. Design Concept**
- 006 II. Surgery
- 007 III. Prosthetics
- 010 IV. Maintenance
- ⁰¹¹ Fixture Product & Packaging
- 011 **I. Dimension**
- 012 II. Fixture Size
- 016 III. Packaging
- 017 Cover Screw & Healing Abutment
- 018 Abutment & Prosthetic options
- 018 I. Fixture level Prosthesis
- nge 1. Extra EZ Post
- 034 II. Abutment level Prosthesis
 - 1. Solid Abutment & Components
- 038 2. Octa Abutment & Components
- 3-1. Multi-unit Abutment & Components
- 059 3-2. Multi-unit Abutment & Components
- of 2 4. Flat Abutment & Components
- 066 III. Overdenture Prosthesis
 - 1. MegaGen Overdenture System
- 067 2. Meg-Loc Abutment & Component
- o70 3. Meg-Ball Abutment & Component
- 073 4. Meg-Magnet Abutment & Component
- 5. Meg-Rhein Abutment & Components

- 082 AnyRidge Kit
- 082 I. Abutment Selection Guide Kit
- 083 II. Surgical Kit
- 089 III. Prosthetic Kit
- 091 IV. Bone Profiler Kit
- 092 V. Optional components
 - R2GATE Full Surgica
- 094 I. R2GATE Full Surgical Kit
- 100 R2GATE Universal K
- 106 Digital Material
- 106 I. ZrGEN
- 106 II. TIGEN
- 121 Clinical Cases

Product coordinator: Hyun wook An, rnd_dir@imegagen.com

Characteristics & Advantages

I. Design Concept

Small but Strong Abutment Screw

Diameter 1.8mm

Thicker abutment wall

More favorable for path adjustment

Optimum post taper

- ✓ Different taper according to post diameter (8°, 10°, 12°, 14°)
- ✓ Wider diameter has more taper!

Various cuff height 2.0, 3.0, 4.0, 5.0mm

Single connection

Abutment can be used on any size of fixture

Wider fixture in a narrow crest Maximizes long-term survival of a fixture

No cutting edge, but strong self-threading

- ✓ Sharp cutting flutes slice and widen bone gradually.
- ✓ No wobbling on cortical slope in anterior immediate placement case.

Narrow diameter

Wider fixture into a small osteotomy socket (less invasive surgery) is important to preserve marginal hard and soft tissues.

Various post height 4.0, 5.5, 7.0, 9.0mm

Flexibility of Imm

1mm trimmable margin gives restorative flexibility without changing B-L dimension

Biologic S-line

Biologic S-line provides seamless natural-looking emergence profile

Maximum preservation of cortical bone

Important for esthetics and long-term prognosis

Knife-Threads

Offer progressive bone condensing, ridge expansion, maximized compressive force resistance, and minimized shear force production.

Taper design

Easy to place always guarantees excellent initial stability

AnyRidge Fixture Line up

Same core diameter, but different thread depth

Core	Fixture Diameter									
Diameter	Ø 3. 5	Ø 4.0	Ø 4. 5	Ø 5.0	Ø 5. 5	Ø 6.0	Ø 6. 5	Ø 7.0	Ø 7. 5	Ø8.0
∅2.8										
Thread depth	0.3									
∅3.3										
Thread depth		0.35	0.6	0.85	1.1					
∅3.8										
Thread depth			0.35	0.6	0.85					
∅4.0										
Thread depth				0.45	0.7	0.95				
∅4.3										
Thread depth				0.35	0.6	0.85				
∅4.8									-	
Thread depth					0.35	0.6	0.85	1.1	1.35	1.6

II. Surgery

Excellent initial stability, even at compromised bone density.

AnyRidge® Fixture cuts bone smoothly and condenses it simultaneously.

1. Fixture placement


Soft bone

The super self-tapping threads have a single core diameter that facilitates minimal site preparation by utilizing a smaller osteotomy to place a wider fixture with special threads.

Hard bone

AnyRidge* Fixture with its super self-tapping thread design is easier than other traditional implants at hard bone.

*Caution! : The osteotomy socket (drilling) size should almost reach the size of fixture to avoid getting stuck in the bone during placement.



Easy way to avoid stuck in the bone during AnyRidge implant placement

1. Due to extremely strong initial stability of AnyRidge fixture, it can be stuck in the middle during placement especially in mandibular hard bone. Please consider 'One millimeter Rule' to avoid this in the best and easiest way. Clinician can customize the drilling sequence once he fully understand the concept and characteristics of AnyRidge system to get preferred initial stability. 'One millimeter Rule' is simple; if an implant engine (40Ncm) stops leaving one millimeter above the crest, use ratchet wrench to screw it down to preferred position. We recommended to place implant platform 0.5~1.0mm under the crest.

2. If a fixture stucks in the middle leaving more than 1mm above the crest in hard mandibular bone, it is recommended to remove it using a wrench rather than trying to screw it down with excessive torque. Please use a cortical bone drill that is included in a surgical kit, the depth of cortical bone drilling can be adjusted according to the bone condition. Then, place the same fixture into the osteotomy socket.

2. Customized drilling Sequence

AnyRidge® system has no fixed drilling protocol, just make your own protocol based on
patient's bone quality to attain preferred initial stability or simply drill an osteotomy
socket to given conditions and then decide the diameter of a fixture.



 Improved drill design has simplified drilling sequence, you can even harvest autogenous bone using these specially designed drills.
 (Recommended speed: 50 RPM, 50 Ncm with saline solution irrigation)

• The best way to get ideal initial stability with AnyRidge system is placing a fixture using a surgical engine, leaving one or two treads above the crest; then use ratchet wrench to place the platform at the desired position.


III. Prosthetics

Better esthetic outcomes from wide variety of prosthetic options! Stop worrying about screw loosening!

 No screw loosening, less biologic width!

• Magic Five (5° Internal connection)

Now you can be free from screw loosening with our unique connection (5 degree morse taper) which gives perfect hermetic sealing. Biologic width is minimized due to no micro gap, and crestal bone health is well maintained.



Performed Retention Test to evaluate the fixture-abutment retention force using Universal Testing Machine -R&D center in MegaGen Implant Co.,Ltd.(2009)-

2. Biologic S-line


Helps to achieve beautiful, natural-looking esthetics.

3. Optimum hex height

Feel AnyRidge connection. It starts with impression taking and lasts until final restoration.


4. All indications, various abutment options

Every case, every shape, every size. Everything was considered to satisfy every need.



>> Two different connections between a fixture & a component

 All transitional and temporary components have 'Ledges' on the bottom



- Cover Screws, Healing Abutments, Impression Coping (transfer and pick-up type), Temporary Cylinders have ledges on the bottom which prevent from cold welding with a fixture.
- Hand Drivers(1.2 Hex) or Impression Drivers can be used easily to screw these components in and out.


2. All permanent abutments will make strong connections with fixtures, even with finger force!



- 25~35Ncm is recommended to connect a permanent abutment into a fixture.
- A fixed abutment cannot be removed with finger force even after complete removal of Abutment Screw because of perfect cold welding. To remove a permanent abutment, Abutment Removal Driver should be used.

How to remove Permanent Abutment from Fixture?

Due to extremely strong connection force, you don't have to worry about screw loosening. Please use our special 'removal driver' when you want to remove an abutment from a fixture.



- 1. Use a Hand Driver(1.2 Hex) to unscrew Abutment Screw.
- 2. Continue to turn counter-clockwisely until you feel a click of disengagement.
- 3. Push down Hand Driver once again to catch and fix Abutment Screw.
- 4. Lift up Hand Driver lightly and continue to turn counter-clockwisely until Abutment Screw engages with the inner screw of Abutment.
- 5. Remove Abutment Screw completely from the abutment
- 6. Insert an 'Abutment Removal Driver' and continue to turn clockwise until the abutment comes out of fixture. You can feel some resistance during screw-down of the Abutment Removal Driver, but don't worry, simple exert is needed disconnect the abutment from the fixture.




IV. Maintenance

Unique and sturdy design provides long-term stability!

1. Higher cortical bone preservation is guaranteed




AnyRidge does not depend on cortical bone for initial stability; decreased stress on cortical bone helps to prevent bone resorption after implantation.



- More cortical bone
- = More soft tissue volume
- = Beautiful gingival line

Advanced coronal design allows maximum cortical bone preservation around implants. Beyond osseointegration, AnyRidge can assure beautiful gingival line by preserving and maintaining more cortical bone.



• A Human Biopsy (2.5 yrs after placement)

The sharp and high alveolar crest (yellow arrow) could be maintained due to biology driven implant design. With this maintenance of alveolar bone, the peri-implant marginal gingiva showed almost no recession during 2.5 years of follow-up, even in the case of limited ridge width.




Fixture Product & Packaging

I. Dimension



Core (mm)	Platform (mm)	Bevel (mm)
Ø3.3	3.5	3.8
<i>V</i> 3.3	3.3	4.0
Ø3.8	4.0	4.5
Ø4.0	4.25	4.75
Ø4.3	4.5	5.0
Ø4.8	5.0	5.5



II. Fixture Size

Small Ø3.5

- Cover Screw included.
- Availability of 7mm product is subject to local approval.

Fixture Diameter (mm)	Core (mm)	Length (mm)	Ref.C
3.5	2.8	7	FANIHX3507C
		8.5	FANIHX3508C
		10	FANIHX3510C
		11.5	FANIHX3511C
		13	FANIHX3513C
		15	FANIHX3515C



Regular Ø4.0

- Cover Screw included.
 Availability of 7mm product is subject to local approval.

Fixture Diameter (mm)	Core (mm)	Length (mm)	Ref.C
	3.3	7	FANIHX4007C
		8.5	FANIHX4008C
4.0		10	FANIHX4010C
4.0		11.5	FANIHX4011C
		13	FANIHX4013C
		15	FANIHX4015C



Regular Ø4.5

- Cover Screw included.
- Availability of 7mm product is subject to local approval.

Fixture Diameter (mm	Core n) (mm)	Length (mm)	Ref.C
		7	FANIHX4507C
		8.5	FANIHX4508C
	3.3	10	FANIHX4510C
	3.3	11.5	FANIHX4511C
		13	FANIHX4513C
4.5		15	FANIHX4515C
4.5	3.8	7	AR384507C
		8.5	AR384508C
		10	AR384510C
		11.5	AR384511C
		13	AR384513C
		15	AR384515C



○ Fixture Size (Continued)

Wide Ø5.0

Cover Screw included



Fixture Diameter (mm)	Core (mm)	Length (mm)	Ref.C
		7	FANIHX5007C
		8.5	FANIHX5008C
	3.3	10	FANIHX5010C
	3.3	11.5	FANIHX5011C
		13	FANIHX5013C
F.O.		15	FANIHX5015C
5.0	3.8	7	AR385007C
		8.5	AR385008C
		10	AR385010C
		11.5	AR385011C
		13	AR385013C
		15	AR385015C

Fixture Diameter (mm)	Core (mm)	Length (mm)	Ref.C
		7	FANIHX5007SC
		8.5	FANIHX5008SC
	4.0	10	FANIHX5010SC
	4.0	11.5	FANIHX5011SC
		13	FANIHX5013SC
5 0		15	FANIHX5015SC
5.0	4.3	7	AR435007C
		8.5	AR435008C
		10	AR435010C
		11.5	AR435011C
		13	AR435013C
		15	AR435015C

Wide Ø5.5

- Cover Screw included



Fixture Diameter (mm)	Core (mm)	Length (mm)	
		7	FANIHX5507C
		8.5	FANIHX5508C
	0.0	10	FANIHX5510C
	3.3	11.5	FANIHX5511C
		13	FANIHX5513C
		15	FANIHX5515C
		7	AR385507C
	3.8	8.5	AR385508C
5.5		10	AR385510C
5.5		11.5	AR385511C
		13	AR385513C
		15	AR385515C
		7	FANIHX5507SC
		8.5	FANIHX5508SC
	4.0	10	FANIHX5510SC
		11.5	FANIHX5511SC
		13	FANIHX5513SC
		15	FANIHX5515SC

Fixture Diameter (mm)	Core (mm)	Length (mm)	
		7	AR435507C
		8.5	AR435508C
	4.3	10	AR435510C
	4.3	11.5	AR435511C
		13	AR435513C
5.5		15	AR435515C
5.5	4.8	7	AR485507C
		8.5	AR485508C
		10	AR485510C
		11.5	AR485511C
		13	AR485513C
		15	AR485515C

Super Wide Ø6.0 - Cover Screw included.

Fixture Diameter (mm)	Core (mm)	Length (mm)	Ref.C
		7	AR406007C
		8.5	AR406008C
	4.0	10	AR406010C
		11.5	AR406011C
		13	AR406013C
	4.3	7	AR436007C
		8.5	AR436008C
6.0		10	AR436010C
		11.5	AR436011C
		13	AR436013C
		7	FALIHX6007C
		8.5	FALIHX6008C
	4.8	10	FALIHX6010C
		11.5	FALIHX6011C
		13	FALIHX6013C



Super Wide Ø6.5 - Cover Screw included.

Fixture Diameter (mm)	Core (mm)	Length (mm)	Ref.C
6.5	4.8	7	FALIHX6507C
		8.5	FALIHX6508C
		10	FALIHX6510C
		11.5	FALIHX6511C
			13



Super Wide Ø7.0 - Cover Screw included.

Fixture Diameter (mm)	Core (mm)	Length (mm)	Ref.C
7.0	4.8	7	FALIHX7007C
		8.5	FALIHX7008C
		10	FALIHX7010C
		11.5	FALIHX7011C
		13	FALIHX7013C



Sixture Size

Super Wide Ø7.5

- Cover Screw included.


Fixture Diameter (mm)	Core (mm)	Length (mm)	Ref.C
	4.8	7	FALIHX7507C
		8.5	FALIHX7508C
7.5		10	FALIHX7510C
		11.5	FALIHX7511C
		13	FALIHX7513C



Super Wide Ø8.0

- Cover Screw included.

Fixture Diameter (mm)	Core (mm)	Length (mm)	Ref.C
	4.8	7	FALIHX8007C
		8.5	FALIHX8008C
8.0		10	FALIHX8010C
		11.5	FALIHX8011C
		13	FALIHX8013C





Fixture Diameter (mm)	Core (mm)	Length (mm)	Ref.C
Ø3.5	3.3		AR333505C
Ø4.0	3.3		AR334005C
CA F	3.3		AR334505C
Ø4.5	3.8		AR384505C
	3.3		AR335005C
Ø5.0	3.8	7	AR385005C
	4.3		AR435005C
	3.3		AR335505C
05.5	3.8		AR385505C
Ø5.5	4.3		AR435505C
	4.8		AR485505C



"Special 7mm" essential for special case Simin Final Fi


For Irregular Ridge


This 'Special 7mm' fixture can be used for non-uniform bone loss case with limited available vertical dimension.

Ø3.5, Ø4.0, Ø4.5, Ø5.0, Ø5.5




III. Packaging






Peel off cover & remove ampule



Separate top¹ & bottom², as shown, to reveal inner ampule with fixture









Connect handpiece to fixture



Make sure fixture is fully connected, then remove from ampule



Place fixture according to drilling sequence



Separate fixture ampule from bottom, as shown, to reveal cover screw holder[®]



Use hand driver to pick up cover screw



Tighten cover screw to fixture

MegaGen ampule! Re-usable as building block *after cleaning and sterilization! less plastic waste!

Cover Screw & Healing Abutment

Cover Screw

- * Included in the fixture package.
- Use with a Hand Driver(1.2 Hex).
- Used for submerged type surgery.
- Protects the inner structure of a fixture.
- Different heights can be chosen according to the position of fixture below the crest.
- 1.6mm and 2.6mm height of Cover Screw can be purchased separately.
- Recommend torque : by hand (5 8Ncm)

Profile Diameter	Height (mm)	Ref.C
	0.8	AANCSF3508
Ø3.5	1.6	AANCSF3516
	2.6	AANCSF3526
Ø5.0	0.5	AANCSF5005
Ø6.0	0.5	AANCSF6005



Healing Abutment

- Use with a Hand Driver(1.2 Hex).
- Used for non-submerged type surgery or for two stage surgery.
- Choose appropriate diameter and height of Healing Abutment according to situation.
- Helps to form suitable emergence profile during period of gingival healing.
- Recommend torque : by hand (5 8Ncm)




Profile Diameter	Height (mm)	Ref.C
	3	AANHAF0403
	4	AANHAF0404
	5	AANHAF0405
Ø4.0	6	AANHAF0406
	7	AANHAF0407
	8	AANHAF0408
	9	AANHAF0409
	3	AANHAF0503
	4	AANHAF0504
	5	AANHAF0505
Ø5.0	6	AANHAF0506
	7	AANHAF0507
	8	AANHAF0508
	9	AANHAF0509
	3	AANHAF0603
	4	AANHAF0604
	5	AANHAF0605
Ø6.0	6	AANHAF0606
	7	AANHAF0607
	8	AANHAF0608
	9	AANHAF0609

Profile Diameter	Height (mm)	Ref.C		
	3	AANHAF0703		
	4	AANHAF0704		
	5	AANHAF0705		
Ø7.0	6	AANHAF0706		
	7	AANHAF0707		
	8	AANHAF0708		
	9	AANHAF0709		
	3	AANHAF0803		
	4	AANHAF0804		
	5	AANHAF0805		
Ø8.0	6	AANHAF0806		
	7	AANHAF0807		
	8	AANHAF0808		
	9	AANHAF0809		
	3	AANHAF1003		
	4	AANHAF1004		
	5	AANHAF1005		
Ø10.0	6	AANHAF1006		
	7	AANHAF1007		
	8	AANHAF1008		
	9	AANHAF1009		

Abutment & Prosthetic Options

I. Fixture Level Prosthesis




Impression Copings

Impression Coping

(2-piece, Transfer Type) (For Closed-tray Technique)

- Streamlined shape ; easy to transfer.
- Anti-rotation grooves match with hex structure of fixtures.
- Should be tightened with Impression Driver (Page.352)
- Special impression coping screw which can be used with a 1.2mm hex driver is available on request.

Profile Diameter	Height (mm)	Туре	Ref.C
~	12		AANITH4012T
Ø4.0	16	O Diago	AANITH4016T
Ø5.0	12	2-Piece	AANITH5012T
Ø5.0	16		AANITH5016T
04.0	12	2-Piece Hand driver (1.2 Hex)	AANITH4012HT
Ø4.0	16		AANITH4016HT
ØF 0	12		AANITH5012HT
Ø5.0	16	, í	AANITH5016HT




Impression Coping

(1-piece, Transfer Type) (For Closed-tray Technique)

- Should be tightened with Impression Driver (Page.352)
- Special impression coping screw which can be used with a 1.2mm hex driver is available on request.

	Profile Diameter	Height (mm)	Туре	Ref.C
	Ø4.0	12		AANITN4012
	04.0	16	4 D:	AANITN4016
	ØF 0	12	1-Piece	AANITN5012
	Ø5.0	16		AANITN5016
	~	12	1-Piece Hand driver (1.2 Hex)	AANITN4012H
	Ø4.0	16		AANITN4016H
	Ø5.0	12		AANITN5012H
		16		AANITN5016H




Impression Coping

(2-piece, Pick-up Type) (For Open-tray Technique)

- Guide Pins : AANGPP0010 (7mm : Short) / AANGPP0015 (12mm : Long) / AANGPP0020 (20mm : Extra-long)

- Square structure ; strong anti rotation function.
- Designed for easy and accurate pick-up impression.
- Extra-long guide pin can be purchased separately.

Profile Diameter	Height (mm)	Туре	Ref.C
	12		AANIPH4012T
Ø4 0	16		AANIPH4016T
<i>1</i> 04.0	12	2-Piece	AANIPN4012T
	16		AANIPN4016T
	7		AANIPH5007T
QE O	12		AANIPH5012T
Ø5.0	7		AANIPN5007T
	12		AANIPN5012T



Lab Analog & Temporary Abutments

Lab Analog

Profile Diameter	Color	Ref.C
Ø3.5	Magenta	AANLAF35
Ø4.0 ~ Ø5.5	Blue	AANLAF4055
Ø6.0 ~ Ø8.0	Yellow	AALLAF6080




Temporary Abutment

(Titanium)

- Multi Post Screw(AANMSF) included.
- For making provisional restoration.
- Grooved on the post allows strong resin adherence.
- Recommend torque: 25Ncm

Profile Diameter	Cuff Height (mm)	Туре	Ref.C
Ø4.0	2	Hex	AANTMH4012T
	2	Non-Hex	AANTMN4012T




Fuse Abutment


- Straight, 15°, 25°; Multi Post Screw(AANMSF) included + Fuse Cap included.
- Milling; Multi Post Screw(AANMSF) included.
- Recommend torque : 25Ncm



NEW: 4mm cuff height available

→ Adequate for deeply placed implants or thick gingival cases






→ Fuse Abutment




Why is the 'Fuse Abutment' essential partner for a temporary crown?

Design concept of Fuse Abutment™






Rationale of Fuse Abutment™



Performed compressive strength test to evaluate the micro movement for bone density using universal testing machine
-R&D center in Meaagen Implant Co..Ltd.(2012)-

Compressive strength test of Fuse Abutment



Performed compressive strength test to evaluate the yield strength for Fuse Abutment using universal testing machine -R&D center in Megagen Implant Co. Ltd (2012)- In 1992, Brunski JB. reported that the implant may has a higher possibility of fibrointergration than osseointegration between bone and implant surface when movements of more than 100 um occur on the fixture during osseointegration period. (John B. Brunski, Biomechanical factors affecting the bone-dental implant interface. Clinical Materials, Vol. 10, 153-201) Therefore, the implant was needed to protected not to move when immediate loading is carried out. However, it is not easy to manage loading on the fixture, even when we used a resin temporarily with a titanium cylinder. It was thought that it was partly because of the metal component of temporary cylinder, which can deliver excessive forces to the fixture. This was one of the reasons which made clinicians hesitate the immediate loading procedure. So it was necessary to develop a special temporary cylinder. It should have been broken under the force which could lead fibrointegration or failure of osseointegration to protect the fixture. and it would be preferred if it was easy to make a temporary crown on this particular temporary cylinder. We tried to measure the force causing movement

of 100µm on a fixture which was placed securely into adequate density of bone without defect. First, AnyRidge implants were placed into the internationally recognized standard bone block with more 40Ncm torque force and an abutment was connected on each implant. Instron equipment was used to measure the force to move a fixture 100µm. The average force was 220N (22.4 kgf). Therefore, if the new temporary abutment can be fractur under this force, it might protect the fixture from movement or failure.




From this experiment, we could developed a special temporary abutment which has lower fracture threshold of less than 200 N (20.4 kgf). It was named as Fuse Abutment. Also it has an anatomic profiles to make temporary prosthetics more esthetic.

○ Abutment Options (Continued)

EZ Post Abutment

- Multi Post Screw(AANMSF/AANMST) included.
- Use with a Hand Driver (1.2 Hex).
- · Esthetic gold coloring.
- Two different post heights. (5.5, 7.0mm)
- Four different profile diameters. (Ø4.0, 5.0, 6.0, 7.0)
- Four different cuff heights. (2.0, 3.0, 4.0, 5.0mm)
- Recommend torque: 35Ncm



Profile Diameter	Cuff Height(mm)	Post Height(mm)	Туре	Ref.C
	2			AANEPH4025L
	3			AANEPH4035L
	4	5.5		AANEPH4045L
04.0	5		Lieu	AANEPH4055L
Ø4.0	2		Hex	AANEPH4027L
	3	7		AANEPH4037L
	4	7		AANEPH4047L
	5			AANEPH4057L
	2			AANEPN4025L
	3	F F		AANEPN4035L
	4	5.5		AANEPN4045L
04.0	5		Niere Lieur	AANEPN4055L
Ø4.0	2		Non-Hex	AANEPN4027L
	3	7		AANEPN4037L
	4	7		AANEPN4047L
	5			AANEPN4057L
	2			AANEPH5025L
	3	<i></i>		AANEPH5035L
	4	5.5		AANEPH5045L
Ø5.0	5		Hex	AANEPH5055L
25.0	2		пех	AANEPH5027L
	3	7		AANEPH5037L
	4	1		AANEPH5047L
	5			AANEPH5057L
	2			AANEPN5025L
	3	5.5		AANEPN5035L
	4	0.0		AANEPN5045L
Ø5.0	5		Non-Hex	AANEPN5055L
Ø3.0	2		NOH-HEX	AANEPN5027L
	3	7		AANEPN5037L
	4			AANEPN5047L
	5			AANEPN5057L


Profile Diameter	Cuff Height(mm)	Post Height(mm)	Туре	Ref.C
	2			AANEPH6025L
	3	F F		AANEPH6035L
	4	5.5		AANEPH6045L
Ø6.0	5		Hov	AANEPH6055L
Ø6.0	2		Hex	AANEPH6027L
	3	7		AANEPH6037L
	4	7		AANEPH6047L
	5			AANEPH6057L
	2			AANEPN6025L
	3			AANEPN6035L
	4	5.5		AANEPN6045L
6 0.0	5			AANEPN6055L
Ø6.0	2		Non-Hex	AANEPN6027L
	3	7		AANEPN6037L
	4	7		AANEPN6047L
	5			AANEPN6057L
	2			AANEPH7025L
	3			AANEPH7035L
	4	5.5		AANEPH7045L
07.0	5			AANEPH7055L
Ø7.0	2		Hex	AANEPH7027L
	3	7		AANEPH7037L
	4	7		AANEPH7047L
	5			AANEPH7057L
	2			AANEPN7025L
	3			AANEPN7035L
	4	5.5		AANEPN7045L
07.0	5			AANEPN7055L
Ø7.0	2		Non-Hex	AANEPN7027L
	3	-		AANEPN7037L
	4	7		AANEPN7047L
	5			AANEPN7057L

○ Abutment Options (Continued)

Gold Abutment

- Multi Post Screw(AANMSF/AANMST) included.
- · Useful to make a customized abutment in difficult situations.


 Precious and non-precious alloys.
- Melting point of gold alloy: 1063°C
- Threaded sleeves for convenient Resin / Wax-up.
- Recommend torque : 30Ncm



Zirconia Abutment

- Multi Post Screw(AANMSF/AANMST) included.
- · For esthetic use.
- · Natural white color with pre-sintered zirconia sleeve.
- · Presinpered Zircornia Abutment.
- Preparable at the chair side with a diamond bur.
- Recommend torque: 35Ncm


Profile Diameter	Cuff Height(mm)	Post Height(mm)	Туре	Ref.C
Ø4.0		11	L.Levi	AANZAH4012L
Ø5.0		11	Hex	AANZAH5012L



Milling Abutment

- Multi Post Screw(AANMSF/AANMST) included.
- · Long post enables easier customization from milling.
- Recommend torque: 35Ncm


	Profile Diameter	Cuff Height(mm)	Post Height(mm)	Ref.C
		2		AANMAH4029L
	Ø4.0	3	9	AANMAH4039L
		4	9	AANMAH4049L
		5		AANMAH4059L
		2		AANMAH5029L
	ØF 0	3	0	AANMAH5039L
	Ø5.0	4	9	AANMAH5049L
		5		AANMAH5059L
		2		AANMAH6029L
	Ø0.0	3	0	AANMAH6039L
	Ø6.0	4	9	AANMAH6049L
	07.0	5		AANMAH6059L
		2		AANMAH7029L
		3	0	AANMAH7039L
	Ø7.0	4	9	AANMAH7049L
		5		AANMAH7059L



Milling Abutment Type II (BOPT Abutment)

- AnyRidge Internal : Multi Post Screw (AANMSF/ AANMST) included.
- Long post enables easier customization from milling.
- Recommend torque: 35Ncm

Profile Diameter	Cuff Height(mm)	Post Height(mm)	Туре	Ref.C
04.0		40	Hex	ARBOT4019HL
Ø4.0	I	19	Non-Hex	ARBOT4019NL



B.O.P.T (Biologically OrientedPreparation Technique)

MegaGen family thanks to MD. Oscar Alonso Gonzalez & Dr. Fabio Galli for the suggestion of B.O.P.T abutment

 To obtain thick, healthy and stable soft tissue annound tooth




Fig 1. With its divergent profile, it tends to stabilize the circular fibers of the connective tissue towards apical.


Fig 2. In the same way as with the teeth, this abutment facilitates the stabilization of the circular fibers of the connective tissue at a more coronal level compared to a standard rehabilitation.

Characteristics of B.O.P.T

- 1. Morphology without a finish line.
- 2. Conical Shape.
- 3. Prosthetic Platform Switching

B.O.P.T Clinical Case

- Courtesy of Dr. Fabio Galli



Implant Placement

Post-op panoramic view


Final Restoration

4 years F/U

○ Abutment Options (Continued)

Angled Abutment

- Multi Post Screw(AANMSF/AANMST) included.
- Two different angulations. (15°, 25°)
- Four different profile diameters. (Ø4.0, 5.0, 6.0, 7.0)
- Four different cuff heights. (2, 3, 4, 5mm)
- Can cover 12 different directions.
 [six to the surface(Hex), six to the edge of hex(Hex-E)]
- Esthetic gold coloring.
- Minimized screw head length needs minimum height to prevent milling problems.
- Recommend torque: 35Ncm





	Profile Diameter	Cuff Height(mm)	Post Height(mm)	Туре	Angle	Ref.C
		2				AANAAH4215L
		3		Hex		AANAAH4315L
		4		пех	15°	AANAAH4415L
		5				AANAAH4515L
		2			15	AANAAE4215L
		3		Hex-E		AANAAE4315L
		4		I IEX-L		AANAAE4415L
	Ø4.0	5	7			AANAAE4515L
	04.0	2	<i>'</i>			AANAAH4225L
		3		Hex		AANAAH4325L
		4		пех		AANAAH4425L
		5			050	AANAAH4525L
		2		Hex-E	25°	AANAAE4225L
		3				AANAAE4325L
		4				AANAAE4425L
		5				AANAAE4525L
		2		Hex		AANAAH5215L
		3				AANAAH5315L
		4			15°	AANAAH5415L
		5				AANAAH5515L
		2			15	AANAAE5215L
		3		Hex-E		AANAAE5315L
		4		⊓ex-⊏		AANAAE5415L
	OF O	5	7			AANAAE5515L
	Ø5.0	2	/			AANAAH5225L
		3		Llov		AANAAH5325L
		4		Hex		AANAAH5425L
		5			25°	AANAAH5525L
		2			25	AANAAE5225L
		3		Hex-E		AANAAE5325L
		4		пех-Е		AANAAE5425L
		5				AANAAE5525L

Profile Diameter	Cuff Height(mm)	Post Height(mm)	Туре	Angle	Ref.C
	2				AANAAH6215L
	3		Hex		AANAAH6315L
	4		пех		AANAAH6415L
	5			15°	AANAAH6515L
	2			15	AANAAE6215L
	3		Hex-E		AANAAE6315L
	4		I IEX-L		AANAAE6415L
Ø6.0	5	7			AANAAE6515L
20.0	2	,			AANAAH6225L
	3		Hex		AANAAH6325L
	4		LICX		AANAAH6425L
	5			25°	AANAAH6525L
	2			20	AANAAE6225L
	3		Hex-E		AANAAE6325L
	4				AANAAE6425L
	5				AANAAE6525L
	2				AANAAH7215L
	3		Hex		AANAAH7315L
	4		LICX		AANAAH7415L
	5			15°	AANAAH7515L
	2			13	AANAAE7215L
	3		Hex-E		AANAAE7315L
	4		I IEX-L		AANAAE7415L
Ø7.0	5	7			AANAAE7515L
Ø1.0	2	,			AANAAH7225L
	3		Hex		AANAAH7325L
	4	LICX		AANAAH7425L	
	5			25°	AANAAH7525L
	2			25	AANAAE7225L
	3		Hex-E		AANAAE7325L
	4		ITEX-E		AANAAE7425L
	5				AANAAE7525L


CCM Abutment

- Multi Post Screw(AANMSF/AANMST) included.
- Useful to make a customized abutment in difficult situations.
- Can be casted with non-precious alloys (Ni-Cr, Cr-Co alloys).
- Non-precious melting temperature : Depend on Manufacturer
- Threaded sleeves for convenient Resin / Wax-up.
- Melting temperature of CCM : 1300~1400 °C
- Recommend torque : 35Ncm

Hex AANCAH4012L Non-Hex AANCAN4012L Platform level	Profile iameter	Cuff Height(mm)	Post Height(mm)	Туре	Ref.C
Non-Hex AANCAN4012L Platform	Ø4.0	1	11	Hex	AANCAH4012L
	<i>0</i> 4.0	' '	11	Non-Hex	AANCAN4012L

I. Fixture Level Prosthesis

1. Fixture Level Prosthesis_Extra EZ Post



S2 Option for successful 'second molar implant'

AnyRidge challenges to the HIGH SURVIVAL RATE even at the second molar

You may already know that

'second molar Implant' has much less success rate than others

1) Simple Literature Reviews:

<General Implant success rate>

99.7% - 10-year survival rate at implant

- van Velzen FJ et al. (2014)

95.6%, 94.4%, 96.1%, 100%, 90.6%, 95.7% - CSR of 759 implants in single-tooth prostheses, cantilever fixed, partial prostheses, fixed partial prostheses, fixed complete prostheses, implant/tooth-connected prostheses, and overdentures - Romeo E et al. (2004)

<Second molar Implant Success rate>

"89.0%" - CSR of 392 implants in the posterior mandible for 6 yrs

- Parein et al. (1997)

"91.1%" - 2" Molar survival Rate for 2 yrs - YK kim et al. (2010)

"82.9%", "91.5%" - Prospective study on 282 implants placed in the Mx and Mn molar position (6 years cumulative study) Becker et al. (1999)

8.16% failure in the Mx, 4.93% in the Mn - Moy et al (2005)

a) Why io

2) Why less success rate at the Second Molar?

Handicaps of the Second Molar Implant;

1. Less quality & quantity of alveolar bone

- Maxillary 2[™] Molar site usually show less quality (Type IV or worse) and/or limited height due to Sinus pneumatization.
- Mandible 2^{et} Molar site usually show less blood supply which is important for adeguate alveolar bone metabolism. And limited height of bone due to the inferior mandibular nerve.

2. Strong Occlusal force

Due to special joint system at TMJ, the Second Molar usually endure strong occlusal force during mastication.

3. Hygiene Problem

Due to remote position, it's very difficult to maintain hygiene at the second Molar, especially a t the distal area, So easy to get peri-implantitis than others.

dition

3) How to overcome

less success rate?

Possible solution

- We need an implant system which can provide excellent initial stability[®] even at the loose bone and limited height of bone.
- We need an implant system which can provide enough surface area[®] for osseointegration, even at the limited height of bone.
- We need to provide enough space for angiogenesis and blood supply[®] for more active bone remodeling.

We need stronger implant fixture and abutment connection[®] to withstand occlusal forces and lateral movement.

We need to choose **adequate material**⁵ for abutment and crown, which retains much less plaque, even with less accessibility and hygiene skills.

4) MegaGen's suggestion for the second molar implant


"S2 Option" strongly recommended by KOLs of MegaGen.


① Excellent initial stability at loose bone


2 Enough surface area for osseointegration

- Already well-know advantages of AnyRidge Implant System.









[ISQ value comparison between AnyRidge and EZ plus]

You may already know that

'second molar Implant' has much less success rate than others

2" molar immediate placement case

Fig.1 After extraction

Fig.2 Immediate placement

Fig.3 Final restoration

Fig.4 7yr F/U


[EZ Plus]

- 4.8 core diameter and deep thread AnyRidge implant will create very strong and satisfactory initial stability at the large extraction socket of second molar.

- Courtesy of Dr. Kwang Bum Park

3 Enough space for angiogenesis and blood supply through the inter-thread space

- Knife thread design of AnyRidge implant creates the maximum space for blood supply



[AnyRidge]

4 Stronger fixture and abutment connection


Fixture Selection

- For the strength of successful second molar implant, our KOLs strongly recommend to use 'Core Diameter' widen than 3.8mm.
- If there is enough width of bone, 4.3mm or 4.8mm core AnyRidge fixture would be better.
- At the large extraction socket of second molar, we recommend 4.8mm core and deep thread AnyRidge implant.

Core					Fixture	Diameter				
Diameter	Ø3.5	Ø4.0	Ø4.5	Ø5.0	Ø5.5	Ø6.0	Ø6.5	Ø7.0	Ø7.5	Ø8.0
Ø 2.8										
Thread depth	0.3									
Ø 3.3										
Thread depth		0.35	0.6	0.85	1.1					
Ø 3.8										
Thread depth			0.35	0.6	0.85					\$2 Option
∅4.0										
Thread depth				0.45	0.7	0.95				
Ø 4.3										
Thread depth				0.35	0.6	0.85				
Ø 4.8							-		-	#
Thread depth					0.35	0.6	0.85	1.1	1.35	1.6

Abutment Selection

- 5° AnyRidge connection is really strong and shows almost no biological width.
- Double offset (Implant switching and Abutment switching) is very helpful to improve soft tissue esthetics and health.
- However, at the second molar implant, the strength against lateral occlusal force is more critical than esthetics.
- So our KOLs strongly recommend to use 'Extra EZ Connection' for abutment.



This 'Double connection' has double advantages.

- 1. Strong resistance to lateral occlusal forces
- 2. No sinking of prosthetics
 - Most of internal connection shows 30~50 µm of sinking following delivery of crown
 - S2 Option will not show sinking phenomenon, while maintaining the 5° internal connection

⑤ Adequate material for hygiene

Our KOLs recommend zirconia customized abutment and/or zirconia monolithic crown for the second molar implant.



ZrGEN(Extra EZ)


ZrGEN is the brand name of Mega-Gen Titanium Base. The strength of ZrGEN frees you from the chipping to conventional PFM prosthesis. Monolithic zirconia crowns have no metal substructure,

enhancing better survival rate!


Bacterial Adhesion on Commercially Pure Titanium and Zirconium Oxide Disks: An In Vivo Human Study
Antonio Scarano, Maurizio Piattelli, Sergio Caputi, Gian Antonio Favero, and Adriano Piattelli JP 2004

The mucosal barrier at implant abutments of different materials

Maria Welander, Ingemar Abrahamsson, Tord Berglundh COIR19, 2008; 635-641



Titanium. A homogeneous layer of cocci or filamentous bacteria covers the titanium surface




Zirconium oxide. A small number of bacteria cover the zirconium oxide surface.



(from left: Ti, ZrO2, Ti, Au/Pt-alloy) in place 1 month after implant placement

- However, the Zirconia customized abutment has limitations on strength which leads fracture of zirconia and/or cement-break between ti-insert and Zirconia abutment.
- So MegaGen developed new 'Ti-insert' for the stronger customized abutment!




S2 Option with AnyRidge Clinical Case

Clinical Case 1

- Courtesy of Dr. Seung Yeup Lee

S2 Option Line-up with AnyRidge implant can be the best solution in posterior zone

- Fig 1. Initial Photo
- Fig 2. Harvest Autogenous Bone
- Fig 3. Implant placement
- Fig 4. Provisionalization
- Fig 5, 6. Before / After Surgery
- Fig 7. Final Delivery
- Fig 8. 6 yrs F/U



Clinical Case 2

- Courtesy of Dr. Seuna Yeup Lee

S2 Option Line-up with AnyRidge implant can be the best solution in posterior zone


- Fig 1. Intra Oral before surgery
- Fig 2. Panorama view
- Fig 3. After Implant Placement
- Fig 4. Connect Extra EZ-Post
- Fig 5. Zirconia Customized Abutment using ZrGen
- Fig 6, 7. Connect PMMA



Extra EZ Post Abutment

Extra EZ Post Abutment

- Multi Post Screw(AANMSF/AANMST) included.
- Useful when fixture is exposed over the gum line.
- Recommend torque : 35Ncm



EZ Post Type

	or iyb	_				
Core Diameter	Profile Diameter	Cuff	Ty	ype	Ref.C	
		2			ARNEEH5025L	
		3		Hex	ARNEEH5035L	
		4		пех	ARNEEH5045L	
	Ø5.0	5			ARNEEH5055L	
	<i>1</i> 05.0	2			ARNEEN5025L	
		3			ARNEEN5035L	
		4		Non-Hex	ARNEEN5045L	
Ø3.3		5	Bevel		ARNEEN5055L	
<i>\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\</i>		2			ARNEEH6025L	
		3		Hex	ARNEEH6035L	
		4		пех	ARNEEH6045L	
	Ø6.0	5			ARNEEH6055L	
	<i>1</i> 00.0	2			ARNEEN6025L	
		3		Non-Hex	ARNEEN6035L	
		4		Non-Hex	ARNEEN6045L	
		5			ARNEEN6055L	


Core Diameter	Profile Diameter	Cuff	Ту	/ре	Ref.C
		2			ARREEH6025L
		3		Hex	ARREEH6035L
		4		пех	ARREEH6045L
	Ø6.0	5			ARREEH6055L
	00.0	2			ARREEN6025L
		3		Non-Hex	ARREEN6035L
		4	Platform		ARREEN6045L
040		5			ARREEN6055L
Ø4.0		2		Hex	ARREEH7025L
		3			ARREEH7035L
		4		пех	ARREEH7045L
		5			ARREEH7055L
	Ø7.0	2			ARREEN7025L
		3		N	ARREEN7035L
		4		Non-Hex	ARREEN7045L
	-	5			ARREEN7055L

Milling Type

	, ,,				
Core Diameter	Profile Diameter	Cuff Height	Post Height	Type	Ref.C
Ø3.3	Ø5.5			Bevel	AANEEH3335L
Ø4.0	Ø5.5	3	5.5	Dietferre	AANEEH4035L
Ø4.8	Ø6.5			Platform	AANEEH4835L



- AANEEH4035 used for fixture (Ø5.0, Ø5.5_Core ø4)
 AANEEH4035 is for the Core Diameter ø4.0mm (Fixture Diameter Ø5.0~5.5mm). It also can be used for Fixture Diameter Ø6.0~8.0mm for platform switching.
- AANEEH4835 used for fixture (Ø6.0~8.0) Recommend torque : 35Ncm



Components for Extra EZ Post Abutment

Cover Screw

(Extra Type)

- Included in the fixture package.
- Use with a Hand Driver(1.2 Hex).
- · Used for submerged type surgery.
- Protects the inner structure of a fixture.
- Different heights can be chosen according to the position of fixture below the crest.
- 1.6mm and 2.6mm height of Cover Screw can be purchased separately.
- Recommend torque : by hand (5 8Ncm)

Core Diameter	Profile Diameter	Туре	Ref.C
Ø3.3	Ø4.0	Bevel	AANCSF4008
Ø4.0	Ø4.25	Platform	AANCSF4208



Extra Healing Abutment

- Use with a Hand Driver(1.2 Hex).
- Used for non-submerged type surgery or for two stage surgery.
- Choose appropriate diameter and height of Heal- ing Abutment according to situation.
- Helps to form suitable emergence profile during period of gingival healing.
- Recommend torque : by hand (5 8Ncm)

Ø3.3 ARNEHA503 ARNEHA504 ARNEHA504 ARNEHA505 ARNEHA506 ARNEHA506 ARNEHA507 ARNEHA603 ARNEHA603 ARNEHA604 ARNEHA605 ARNEHA605 ARNEHA606 ARNEHA606 ARREHA403 ARREHA403 ARREHA404 ARREHA406 ARREHA406 ARREHA407 ARREHA603 ARREHA604 ARREHA606 ARREHA605 ARREHA606 ARREHA606 ARREHA700 ARREHA704 ARREHA706 ARREHA706 ARREHA706 ARREHA706 ARREHA706 ARREHA707 AANHAF484	Diameter	Diameter	(mm)	Type	Ref.C
Ø3.3 Bevel ARNEHA505 ARNEHA507 ARNEHA507 ARNEHA603 ARNEHA603 ARNEHA604 ARNEHA605 ARNEHA605 ARNEHA606 ARNEHA606 ARNEHA403 ARREHA403 ARREHA404 ARREHA404 ARREHA405 ARREHA406 ARREHA406 ARREHA606 ARREHA606 ARREHA607 ARREHA606 ARREHA606 ARREHA606 ARREHA606 ARREHA703 ARREHA703 ARREHA706 ARREHA706 ARREHA706 ARREHA706 ARREHA706			3		ARNEHA503
Ø3.3 Bevel ARNEHA506 ARNEHA607 ARNEHA603 ARNEHA604 ARNEHA606 ARNEHA606 ARNEHA606 ARNEHA606 ARNEHA606 ARNEHA400 ARREHA403 ARREHA404 ARREHA405 ARREHA406 ARREHA406 ARREHA606 ARREHA606 ARREHA606 ARREHA606 ARREHA606 ARREHA606 ARREHA606 ARREHA607 ARREHA606 ARREHA703 ARREHA703 ARREHA705 ARREHA706 ARREHA706 ARREHA706 ARREHA706 ARREHA706			4		ARNEHA504
Ø3.3 Bevel ARNEHA507 ARNEHA603 ARNEHA603 ARNEHA604 ARNEHA605 ARNEHA605 ARNEHA606 ARNEHA607 ARREHA403 ARREHA404 ARREHA404 ARREHA405 ARREHA406 ARREHA406 ARREHA406 ARREHA606 ARREHA603 ARREHA603 ARREHA604 ARREHA606 ARREHA606 ARREHA703 ARREHA703 ARREHA704 ARREHA705 ARREHA706 ARREHA706 ARREHA706 ARREHA706		Ø5.0	5		ARNEHA505
Ø3.3 Bevel ARNEHA603 ARNEHA604 ARNEHA604 ARNEHA605 ARNEHA605 ARNEHA606 ARNEHA607 ARREHA403 ARREHA403 ARREHA404 ARREHA405 ARREHA405 ARREHA406 ARREHA406 ARREHA603 ARREHA603 ARREHA603 ARREHA604 ARREHA606 ARREHA606 ARREHA606 ARREHA703 ARREHA703 ARREHA704 ARREHA705 ARREHA706 ARREHA706 ARREHA706 ARREHA706 ARREHA706 ARREHA707			6		ARNEHA506
3 ARNEHA603 ARNEHA604 ARNEHA605 ARNEHA606 7 ARREHA403 ARREHA403 ARREHA404 ARREHA406 ARREHA406 ARREHA406 ARREHA406 ARREHA406 ARREHA407 ARREHA603 ARREHA603 ARREHA604 ARREHA606 ARREHA606 ARREHA606 ARREHA606 ARREHA606 ARREHA606 ARREHA703 ARREHA703 ARREHA706 ARREHA706 ARREHA706 ARREHA706 ARREHA706	Ø0.0		7	Doval	ARNEHA507
Ø6.0 5 6 ARNEHA605 ARNEHA606 ARNEHA606 ARNEHA607 ARREHA403 ARREHA403 ARREHA404 ARREHA405 ARREHA405 ARREHA406 ARREHA603 ARREHA603 ARREHA604 ARREHA606 ARREHA606 ARREHA606 ARREHA703 ARREHA703 ARREHA704 ARREHA705 ARREHA706 ARREHA706 ARREHA706 ARREHA707 ARREHA707	<i>V</i> 3.3		3	Devei	ARNEHA603
6			4		ARNEHA604
Ø4.2 5 6 7 ARREHA403 ARREHA404 ARREHA405 ARREHA405 ARREHA406 ARREHA406 ARREHA603 ARREHA603 ARREHA606 ARREHA606 ARREHA606 ARREHA607 ARREHA703 ARREHA704 ARREHA705 ARREHA706 ARREHA706 ARREHA706 ARREHA707 ARREHA707		Ø6.0	5		ARNEHA605
3			6		ARNEHA606
Ø4.2 5 6 7 3 4 Ø4.0 5 96.0 5 Platform ARREHA603 ARREHA604 ARREHA605 ARREHA606 ARREHA606 ARREHA703 ARREHA703 ARREHA704 ARREHA705 ARREHA706 ARREHA706 ARREHA707 ARREHA707			7		ARNEHA607
Ø4.2 5 6 7 7 3 4 ARREHA407 ARREHA407 ARREHA603 ARREHA604 ARREHA606 ARREHA606 ARREHA606 ARREHA703 ARREHA703 ARREHA704 ARREHA705 ARREHA706 ARREHA706 ARREHA707 ARREHA707			3		ARREHA403
6 ARREHA406 7 ARREHA407 ARREHA603 ARREHA603 ARREHA604 ARREHA606 ARREHA606 ARREHA606 ARREHA607 ARREHA703 ARREHA703 ARREHA704 ARREHA705 ARREHA706 ARREHA706 ARREHA706 ARREHA707			4		ARREHA404
7 ARREHA407 ARREHA603 ARREHA604 ARREHA606 ARREHA606 ARREHA607 ARREHA607 ARREHA703 ARREHA703 ARREHA704 ARREHA705 ARREHA706 ARREHA706 ARREHA707		Ø4.2	5		ARREHA405
Ø4.0 Ø6.0 5 6 7 3 4 Platform ARREHA604 ARREHA605 ARREHA606 ARREHA607 ARREHA703 ARREHA703 ARREHA704 ARREHA705 ARREHA706 ARREHA706 ARREHA707			6		ARREHA406
Ø4.0 Ø6.0 5 Platform ARREHA604 ARREHA605 ARREHA606 ARREHA606 ARREHA703 ARREHA703 ARREHA704 ARREHA705 ARREHA706 ARREHA706 ARREHA707			7		ARREHA407
Ø4.0 Ø6.0 5 Platform ARREHA605 6 7 ARREHA606 ARREHA607 3 4 ARREHA703 ARREHA704 4 ARREHA705 ARREHA706 6 7 ARREHA707			3		ARREHA603
6 Platform ARREHA606 7 ARREHA607 ARREHA703 ARREHA704 ARREHA705 ARREHA706 7 ARREHA707			4		ARREHA604
6 ARREHA606 7 ARREHA607 3 ARREHA703 4 ARREHA704 Ø7.0 5 ARREHA705 6 ARREHA706 7 ARREHA707	Ø4.0	Ø6.0	5	Dlatform	ARREHA605
3 ARREHA703 4 ARREHA704 Ø7.0 5 ARREHA705 6 ARREHA706 7 ARREHA707			6	Fiationiii	ARREHA606
4 ARREHA704 Ø7.0 5 ARREHA705 6 ARREHA706 7 ARREHA707			7		ARREHA607
Ø7.0 5 ARREHA705 6 ARREHA706 7 ARREHA707			3		ARREHA703
6 ARREHA706 7 ARREHA707			4		ARREHA704
7 ARREHA707		Ø7.0	5		ARREHA705
			6		ARREHA706
Ø4.8 Ø6.5 4 AANHAF484			7		ARREHA707
	Ø4.8	Ø6.5	4		AANHAF484



Lab Analog

Profile Diameter	Color	Ref.C
Ø4.0 ~ Ø5.5	Blue	AANLAF4055
Ø6.0 ~ Ø8.0	Yellow	AALLAF6080






Impression Coping

(2-piece, Transfer Type) (For Closed-tray Technique)

- Streamlined shape; easy to transfer.
- · Anti-rotation grooves match with hex structure of fixtures.
- Should be tightened with Impression Driver (Page.352)
- Special impression coping screw which can be used with a 1.2mm hex driver is available on request.

Height (mm)	Туре	Ref.C
12		AANITH4012T
16	O Diago	AANITH4016T
12	z-Piece	AANITH5012T
16		AANITH5016T
12		AANITH4012HT
16	2-Piece	AANITH4016HT
12	driver	AANITH5012HT
16		AANITH5016HT
	(mm) 12 16 12 16 12 16 12 16 12	(mm) Type 12 16 12 16 12 16 12 16 12 2-Piece Hand driver (1.2 Hex)




Impression Coping

(1-piece, Transfer Type) (For Closed-tray Technique)

- Should be tightened with Impression Driver (Page.352)
- Special impression coping screw which can be used with a 1.2mm hex driver is available on request.

Profile Diameter	Height (mm)	Туре	Ref.C
Ø4.0	12		AANITN4012
<i>1</i> 04.0	16	1-Piece	AANITN4016
Ø5.0	12	1-Piece	AANITN5012
Ø5.0	16		AANITN5016
040	12		AANITN4012H
Ø4.0	16	1-Piece Hand	AANITN4016H
ØF 0	12	driver (1,2 Hex)	AANITN5012H
Ø5.0	16		AANITN5016H




Impression Coping

(2-piece, Pick-up Type) (For Open-tray Technique)

- Guide Pins : AANGPP0010 (7mm : Short) / AANGPP0015 (12mm : Long) / AANGPP0020 (20mm : Extra-long)
- Square structure; strong antirotation function.
- · Designed for easy and accurate pick-up impression.
- Extra-long guide pin can be purchased separately.


Profile Diameter	Height (mm)	Туре	Ref.C
	12		AANIPH4012T
Ø4 0	16		AANIPH4016T
04.0	12		AANIPN4012T
	16	0.0:	AANIPN4016T
	7	2-Piece	AANIPH5007T
QE 0	12		AANIPH5012T
Ø5.0	7		AANIPN5007T
	12		AANIPN5012T



Temporary Abutment (Titanium Extra Type)


- fixture package included.
- Use with a Hand Driver(1.2 Hex).
- · Used for submerged type surgery.
- Protects the inner structure of a fixture.
 Different heights can be chosen according.
- Different heights can be chosen according to the position of fixture below the crest.
- 1.6mm and 2.6mm height of Cover Screw can be purchased separately.
- Recommend torque: by hand (5 8Ncm)

D	Core liameter	Profile Diameter		pe Connection	Ref.C
	Ø3.3	Ø4.5	Dovel	Hex	ARNTAH4510T
	<i>\oldsymbol</i>	W4.5	Bevel	Non-Hex	ARNTAN4510T
	040	04.75	DI-#	Hex	ARRTAH4710T
	Ø4.0	Ø4.75	Platform	Non-Hex	ARRTAN4710T



II. Abutment Level Prosthesis


1. Solid Abutment & Components



Solid Abutment Option

Solid Abutment

- Used in cement retained restoration only.
- Solid Abutment should be placed into patient's mouth before taking impression.
- · Onebody (screw + abutment)
- Should be tightened with a Solid Driver and a Torque Wrench: 35Ncm
- Four different profile diameters. (Ø4.0/5.0/6.0/7.0)
- Should be tightened with special Solid Driver.Wider profile has bigger post angulation.
- (4mm 8°, 5mm 10°, 6mm 12°, 7mm 14°)
- Four different cuff heights. (2/3/4/5mm)
- Three different post heights. (4/5.5/7mm)
- Recommend torque : 35Ncm



Profile Diameter	Cuff Height(mm)	Post Height(mm)	Ref.C
	2		AANSAL4024
	3	4	AANSAL4034
	4	4	AANSAL4044
	5		AANSAL4054
	2		AANSAL4025
Ø4.0	3	5.5	AANSAL4035
	4	5.5	AANSAL4045
	5		AANSAL4055
	2		AANSAL4027
	3	7	AANSAL4037
	4	,	AANSAL4047
	5		AANSAL4057
	2		AANSAL5024
	3	4	AANSAL5034
	4	_	AANSAL5044
	5		AANSAL5054
	2		AANSAL5025
Ø5.0	3	5.5	AANSAL5035
20.0	4	5.5	AANSAL5045
	5		AANSAL5055
	2		AANSAL5027
	3	7	AANSAL5037
	4	,	AANSAL5047
	5		AANSAL5057

Profile Diameter	Cuff Height(mm)	Post Height(mm)	Ref.C
	2		AANSAL6024
	3		AANSAL6034
	4	4	AANSAL6044
	5		AANSAL6054
	2		AANSAL6025
Ø6.0	3	5.5	AANSAL6035
	4	5.5	AANSAL6045
	5		AANSAL6055
	2		AANSAL6027
	3	7	AANSAL6037
	4	' [AANSAL6047
	5		AANSAL6057
	2		AANSAL7024
	3	4	AANSAL7034
	4	4	AANSAL7044
	5		AANSAL7054
	2		AANSAL7025
Ø7.0	3	5.5	AANSAL7035
Ø1.0	4	5.5	AANSAL7045
	5		AANSAL7055
	2		AANSAL7027
	3	7	AANSAL7037
	4	/	AANSAL7047
	5		AANSAL7057

Components for Solid Abutment (Continued)

Snap Impression Coping

- For impression on Solid Abutments.
- 3 colors for different post heights.
- · 4 different diameters for profile diameters. (Ø4, 5, 6, 7)
- · Do not use when abutment is trimmed.

Profile Diameter	Ref.C	, , ,		
	AANSIF440		量量	
Ø4.0	AANSIF455			
	AANSIF470			
	AANSIF540		5 1	
Ø5.0	AANSIF555			
	AANSIF570			
	AANSIF640	_		
Ø6.0	AANSIF655		5-3	
	AANSIF670		==	
	AANSIF740			
Ø7.0	AANSIF755			
	AANSIF770			
				7
				P.D

Comfort Cap

- Protects the Solid Abutment and minimizes irritation to tongue and oral mucosa.
- Can be applied under temporary prosthetics.
 Color coded according to post heights.

Profile Diameter	Post Height(mm)	Ref.C
	4	AANCCF440
Ø4.0	5.5	AANCCF455
	7	AANCCF470
	4	AANCCF540
Ø5.0	5.5	AANCCF555
	7	AANCCF570
	4	AANCCF640
Ø6.0	5.5	AANCCF655
	7	AANCCF670
	4	AANCCF740
Ø7.0	5.5	AANCCF755
	7	AANCCF770

Components for Solid Abutment

Lab Analog

• Directly connected to the Snap Impression Coping in the impression to make a stone model.

Profile Diameter	Height(mm)	Ref.C
	4	AANSLF440
Ø4.0	5.5	AANSLF455
	7	AANSLF470
	4	AANSLF540
Ø5.0	5.5	AANSLF555
	7	AANSLF570
	4	AANSLF640
Ø6.0	5.5	AANSLF655
	7	AANSLF670
	4	AANSLF740
Ø7.0	5.5	AANSLF755
	7	AANSLF770



Burn-out Cylinder


- Fits with a Lab Analog(solid level).
- Easy to wax-up and accurate casting.
- White Cylinder for multiple unit.
- Red Cylinder for single crown.

Profile Diameter	Туре	Ref.C
Ø4.0		AANBCB470
Ø5.0	Multiple	AANBCB570
Ø6.0	Multiple	AANBCB670
Ø7.0		AANBCB770
Ø4.0		AANBCS470
Ø5.0		AANBCS570
Ø6.0		AANBCS670
Ø7.0		AANBCS770



II. Abutment Level Prosthesis

2. Octa Abutment & Components




Components for Octa Abutment (Continued)

Octa Abutment

- Used in manufacturing multiple screw-retained prosthetics.
- Recommend torque : 35Ncm

	ricconninc	na torque.	OOIVOIII
•	Maximum	path angle	: 70°

Profile Diameter	Cuff Height (mm)	Ref.C
	1	AANOAF4010
	2	AANOAF4020
Ø4.0	3	AANOAF4030
	4	AANOAF4040
	5	AANOAF4050
	1	AANOAF0010
	2	AANOAF0020
Ø5.0	3	AANOAF0030
	4	AANOAF0040
	5	AANOAF0050
	1	AANOAF6010
	2	AANOAF6020
Ø6.0	3	AANOAF6030
	4	AANOAF6040
	5	AANOAF6050



Healing Cap

- Cylinder Screw(IRCS200) included.
- · Protects Octa Abutment and minimizes irritation to tongue and oral mucosa.

Profile Diameter	Ref.C
Ø4.0	AANOHC4000T
Ø5.0	IHC400T
Ø6.0	AANOHC6000T




Components for Octa Abutment

Impression Coping

(Transfer)

- Guide Pin(AAOTGP10 / AAOTGP12) included.
- Should be tightened with Impression Driver (Page.256)
- Special impression coping screw which can be used with a 1.2mm hex driver is available on request.


	Profile Diameter	Height (mm)	Туре	Ref.C
		7.5	Octa	AAOITO4010T
	Ø4 0	7.5	Non-Octa	AAOITN4010T
	Ø4.0	0.5	Octa	AAOITO4012T
		9.5	Non-Octa	AAOITN4012T
		7.5	Octa	AAOITO5010T
	ØF 0		Non-Octa	AAOITN5010T
	Ø5.0	9.5	Octa	AAOITO5012T
			Non-Octa	AAOITN5012T
		7.5	Octa	AAOITO6010T
	Ø6.0	7.5	Non-Octa	AAOITN6010T
	Ø6.0 9.5	0.5	Octa	AAOITO6012T
		9.5	Non-Octa	AAOITN6012T



Impression Coping (Pick-Up)

- Guide Pin(AAOPGP10 / AAOPGP12) included.

Profile Diameter	Height (mm)	Туре	Ref.C
	10.0	Octa	AAOIPO4010T
Ø4 0		Non-Octa	AAOIPN4010T
<i>1</i> 04.0	12.0	Octa	AAOIPO4012T
	12.0	Non-Octa	AAOIPN4012T
	05.0 12.0	Octa	AAOIPO5010T
ØF 0		Non-Octa	AAOIPN5010T
25.0		Octa	AAOIPO5012T
		Non-Octa	AAOIPN5012T
	10.0	Octa	AAOIPO6010T
9 0.0	10.0	Non-Octa	AAOIPN6010T
Ø6.0	12.0	Octa	AAOIPO6012T
	12.0	Non-Octa	AAOIPN6012T



Lab Analog

Profile Diameter	Ref.C
Ø3.8	AANOLA4000
Ø4.8	IOA300
Ø5.8	AANOLA6000



Temporary Cylinder

- Cylinder Screw(IRCS200) included.

• Recommend torque : 25Ncm

Profile Diameter	Туре	Ref.C
Ø4 0	Octa	AANOTCO4010T
W4.U	Non-Octa	AANOTCN4010T
Ø5.0	Octa	AANOTCO5010T
Ø5.0	Non-Octa	AANOTCN5010T
000	Octa	AANOTCO6010T
Ø6.0	Non-Octa	AANOTCN6010T








EZ Post Cylinder

- Cylinder Screw(IRCS200) included.

• Recommend torque : 35Ncm

Profile Diameter	Post Height(mm)	Туре	Ref.C
	5.5	Octa	AAOECO4005T
Ø4.0	7.0	Ocia	AAOECO4007T
Ø4.0	5.5	Non-Octa	AAOECN4005T
	7.0	Non-Octa	AAOECN4007T
	5.5	Octa	AAOECO5005T
ØF 0	7.0		AAOECO5007T
Ø5.0	Ø5.0 5.5 Non-Octa	Non Ooto	AAOECN5005T
		AAOECN5007T	
	5.5	Octa	AAOECO6005T
Ø6 0	7.0		AAOECO6007T
Ø6.U	5.5	Non-Octa	AAOECN6005T
	7.0	INOII-OCIA	AAOECN6007T




Components for Octa Abutment

Gold Cylinder

- Cylinder Screw(IRCS200) included.
- For customizing abutment for screw retained multi-unit restoration.
- · Available in both octa(red) and non-octa(white).
- Melting point of gold alloy: 1063°C
- Threaded sleeves allow better retention of resin or wax.
- Available in three diameters (Ø4.0, 5.0, 6.0).
- Recommend torque : 30Ncm


Profile Diameter	Туре	Ref.C
~	Octa	AANGCO4000T
Ø4.0	Non-Octa	AANGCN4000T
Ø5.0	Octa	IOGO100T
	Non-Octa	IOGN100T
Ø6.0	Octa	AANGCO6000T
	Non-Octa	AANGCN6000T



CCM Cylinder

- Cylinder Screw(IRCS200) included.
- Threaded sleeves allow a better retention of resin or wax.
- Available in both octa (pink) and non-octa (yellow) and three diameters (Ø4.0, 5.0, 6.0).
- Recommend torque : 35Ncm
- Melting temperature of CCM : 1300~1400°C
- Can be casted with non-precious alloys (Ni-Cr, Cr-Co alloys).


	Profile Diameter	Туре	Ref.C
	Ø4 0	Octa	AANCCO4000T
	Ø4.0	Non-Octa	AANCCN4000T
	Ø5.0	Octa	AANCCO5000T
	Ø5.0	Non-Octa	AANCCN5000T
	ØC 0	Octa	AANCCO6000T
	Ø6.0	Non-Octa	AANCCN6000T



Plastic Cylinder


- Cylinder Screw(IRCS200) included.
- · Economical option.
- Used for customizing abutment a screw retained multi-unit restorations.
- · Available in both octa (red) and non-octa (white)
- Threaded sleeves allow a better retention of resin or wax.
- Recommend torque: 25Ncm

Profile Diameter	Туре	Ref.C
Ø4.0	Octa	AAOTCO4010T
	Non-Octa	AAOTCN4010T
Ø5.0	Octa	IOPH100T
	Non-Octa	IOPN100T
Ø6.0	Octa	AAOTCO6010T
	Non-Octa	AAOTCN6010T



II. Abutment Level Prosthesis


3-1. Multi-unit Abutment & Components (All-on-4) (N_Type)



>> Multi-unit Abutment™

Multi-unit Abutment Design Concept

MegaGen Implant develops the special abutment named as Multi-unit Abutment, which can be the solution for the edentulous patients. With 4 fixtures placed into patient's ridge and a hybrid denture on those four fixtures, a patient can recover his or her dental condition almost completely. In most cases, Multi-unit Abutments work in a set of 2 x straight type abutment for anterior position and 2 x angled type abutment on posterior position.




Features

You could see how Multi-unit Abutment functions and what benefits you could get from Multi-unit Abutment are as the followings:

- 2 fixtures which are slantly implanted on posterior position are osseointegrated with cancellous bone. These fixtures function as dispersing vertical load on alveolar bone.
- Multi-unit Abutment is only 4 fixtures + 4 abutments. It means that dental surgeon has enough places for surgery. Therefore, it will be easy to find and place 4 fixtures into ridge where abundant cancellous bone exists.
- A doctor can use graft bone material if a patient dosen't have enough alveolar bone. However, the slantly placed fixtures can overcome the patient's insufficient bone by getting good holding strength with this angulation.
- In addition, these angulated fixtures can avoid touching important anatomies, such as mandibular nerve and maxillary sinus.
- All on 4 technique is also possible to do guided surgery using R2GATE Guide with a diagnosis from R2GATE.

Multi-unit Abutment N Type

The solution for the edentulous patients



Benefit

- 1. Easy and economical treatment solution for compromised edentulous cases.
- 2. Expensive and time consuming bone graft may not be necessary.
- 3. Multiple angles (0°, 17°, 30°) support different implant insertion paths.
- 4. Universally compatible with other Multiunit systems.

Available implant System

- AnyRidge Internal
- AnyOne Internal
- AnyOne External

Compatibility with others' Multi-unit level prosthetic components

- √ Post Height
- ✓ Post Diameter
- √ Post Angle
- √ Abutment Angle
- ✓ Cuff Height

Multi-unit Abutment

Multi-unit Abutment [AR] - Straight - MUA Straight Carrier (MUASC) included

- Recommend torque : 35Ncm


Cuff Height (mm)	Туре	Ref.C
1.5	1-piece (M1.8)	MUAARN0015C
2.5		MUAARN0025C
3.5		MUAARN0035C
4.5		MUAARN0045C



Multi-unit Angled Abutment [AR] - 17°

- MUA Screw (MUAARS) included
- MUA Angled Carrier (MUAAC) included
- Recommend torque: 25Ncm


Cuff Height (mm)	Туре	Ref.C
2.5		MUAARH1725LC
3.5	Hex	MUAARH1735LC
4.5		MUAARH1745LC
2.5	Non-Hex	MUAARN1725LC
3.5		MUAARN1735LC
4.5		MUAARN1745LC



Multi-unit Angled Abutment [AR] - 30°

- MUA Screw (MUAARS) included
- MUA Angled Carrier (MUAAC) included
- Recommend torque: 25Ncm

Cuff Height (mm)	Туре	Ref.C
3.5	Hex	MUAARH3035LC
4.5		MUAARH3045LC
3.5	Non-Hex	MUAARN3035LC
4.5		MUAARN3045LC




Components for Multi-unit Abutment (Continued)

Impression coping (Pick-up)

- Guide pin (MUAGP) included
- Use to take an impression at the abutment level.Open tray method.

Connection	Ref.C
Non-Hex	MUAICT



Lab Analog

- · Use to duplicate the Multi-unit abutment in the working
- Available to use as a RP Analog for 3D printed working model.


Head form	Ref.C
Multi-unit Abutment(Nobel)	MUALA



Temporary Cylinder


- Cylinder Screw (MUAS) included
- Use for fabricating acrylic provisional restoration.
 Grooves on the post cylinder allow storing resin adhension.
- · Back-up screw is included.
- · Recommend torque: 15Ncm






CCM Cylinder

- Cylinder Screw (MUAS) 2EA included
- · Use for fabricating screw retained prostheses with metal reinforced or bar structured overdentures.
- Available to cast with non-precious dental alloys (Ni-Cr, Cr-Co alloys)
- Melting temperature of CCM base: 1300~1400°C
- · Back-up screw is included.
- Recommend torque : 15Ncm






Components for Multi-unit Abutment

Healing Cap

- Cylinder Screw (MUAS) 2ea included
- The size of healing cap can be selected depending on soft tissue volume or type of restorations.

Туре	Ref.C
Regular	MUAHCL
Wide	MUAHCWL



Healing Cap Set reference code

Order code : Add "P" after the existing reference code

Ex) MUAHCL → MUAHCP






Try-in Abutment

- Cuff height is indicated with laser markings
- Straight17°, 30°
- · Non-hex type

Angle	Cuff Marking	Ref.C
Straight	1.5 / 2.5 / 3.5 / 4.5	MUTIAAR00C
17°	2.5 / 3.5 / 4.5	MUTIAAR17C
30°	3.5 / 4.5	MUTIAAR30C



Try-in Abutment Set reference code

Order code: MUTIAAR00P



- * Available Systems : AnyRidge, AnyOne Internal, AnyOne External
- * Kit contains Straight, 17° and 30° type of Try-in Abutments (1 each)



Multi-unit Driver

- Use to torque straight type Multi-unit Abutments. Use with a torque wrench (ref code: MTW300A)


Hex	Length	Ref.C	
2.0	10	MUD10	



Right Angle Driver

- Use to torque straight type Multi-unit Abutments.
- Use with latch-type handpiece.
 Use with Meg-TORQ (ref code: MEG_TORQ)


Hex	Length	Ref.C	
2.0	10	MURAD10	



Hand Driver

- Use for abutment screw with 1.2 hex hole.
 Use up to 15° divergent.
 It should use under 30Ncm torque.


Hex	Length	Ref.C
1.2	20	MUHD1220



Removal Driver


- Use for abutment screw with 1.2 hex hole.
 Use up to 15° divergent.
- Exclusively for AnyRidge system.
- It should use under 30Ncm torque.


Hex	Length	Ref.C	
1.2	20	MUARD20	



Screw & Abutment Tightening Torque Guide

- Abutment Screw (M1.8 & M2.0): 25Ncm
- Cylinder Screw (M1.4): 15Ncm
- Straight Abutment (M1.8 &M2.0): 35Ncm






Components for Multi-unit Abutment


Surgical Guide

- The distance between the lines is 7mm
- Put center pin after initial drilling at the centric
 of arch

Angle	Marking Length (mm)	Ref.C
30°	7	MUSG70



>> How to use Surgical Guide



→ Multi-unit Abutment Set Contents

Multi-unit Abutment Healing cap type Set reference code

Order code: Add "HP" after the existing reference code


Ex) MUAARH1725LC → MUAARH1725 HP

Multi-unit Abutment CCM type Set reference code

Order code: Add "P" after the existing reference code


Ex) MUAARH1725LC → MUAARH1725 P






>> Starting Package Contents








Surgical Protocol

Conventional Surgery


1. Initial drilling

For placement of center pin after initial drilling in the centric of the arch. The drilling hole should be in lingual area of the arch to ensure the best result.



2. Guide Bending & Position


Bend according to the patient's arch.




7. Tightening the Abutment

Abutment Screw tightening Torque: 25Ncm

After connecting Abutment Screw, remove Carrier from Abutment. For 17° abutment, you need to tighten it by tilting Driver about 5°. Connect Abutment and check the path using Carrier.




Straight Abutment tightening Torque: 35Ncm After removing Carrier, connect Abutment to the Fixture using Right Angle Driver or MUA Driver.




8. Impression

Take an impression with an individual tray. (Open tray method is strongly recommended to avoid any error in the future.)



9. Healing Cap

Cylinder Screw tightening Torque: 15Ncm
Place Healing Cap on top of Multi-unit abutment, and connect Cylinder Screw with the Hand Driver.




3. Drilling

Drill according to the surgical plan.



4. The fixture is implanted

Place implant fixtures according to the surgical plan and do not exceed torque value (60Ncm)




6. Abutment Selection

Select the appropriate set after checking the angulation and cuff height that were measured with the Try-in abutment. Connect the abutment onto the fixture and check the angulation and the cuff height.



5. Try-in Abutment

Using the laser marking on the Try-in abutment, select the appropriate cuff height and angulation of Multi-unit abutments.



10. Suture




Surgical Protocol

Guided Surgery

1. Guide

Place a R2GATE Guide.



2. Narrow Crest Drill

For the cases with narrow ridge or placing a fixture slanted on the lingual side, you can flatten the surface and drill stably without slipping



8. Setting Temporary and Denture

Reline the temporary denture with resin to fill the space around the Temporary

Cylinder.




7. Connect Temporary Cylinder in the front

Connect the Temporary Cylinders in the front. Make sure that holes in the denture are free from any contact with the Temporary Cylinder. Adjust the holes until there is no contact between the denture and the Temporary Cylinder. *If the Temporary Cylinder is fixed using Guide Pin, resin flow into access hole will be prevented.



9. Connect Temporary Cylinder in the back

Connect rest of the Temporary Cylinders in the back, make sure that the holes in the denture are free from any contact with the Temporary Cylinder. Adjust the holes until there is no contact between the denture and the Temporary Cylinder.



10. Setting Temporary and Denture

All temporary cylinders are picked up in the denture with

resin.




3. Drilling

Drill according to the drilling sequence.



4. Fixture Placement




6. Abutment Selection

Select the appropriate set after checking the angulation and cuff height that were measured with the Try-in abutment. Connect the abutment onto the fixture and check the angulation and the cuff height.



5. Try-in Abutment

Using the laser marking on the Try-in abutment, select the appropriate cuff height and angulation of Multi-unit abutments.





11. Temporary Fixation

Remove Denture and fill up the bottom and other non-resin filled parts with resin and completely fix Temporary Cylinder.



12. Tighten the Denture

Cylinder Screw tightening Torque: 15Ncm Set Denture onto Multi-unit Abutment and tighten cylinder



13. Finish


Close Hole using EZ Seal and finalize the surgery.



▶ Multi-unit Abutment S Type

The solution for the edentulous patients

For the design concept and variable of the Multi-unit Abutment




Benefit

- 1. Retrievability means that doctor can change or retrieve the final prosthetics easily.
- 2. Two types of angulation: 17°, 29°. It means that doctor has various options to angle.
- Various cuff heights (1~5mm): Doctor can have flexiblily on the depth of fixture placement.
- MegaGen's Multi-unit Abutment is perfectly compatible with the prosthetic components of Multi-unit Abutment of 3i implant, and Regular Abutment of Mega-Gen's Exfeel External system.

- AnyRidge Internal
- AnyOne Internal

II. Abutment Level Prosthesis

3-2. Multi-unit Abutment & Components (All-on-4) (S-Type)




Multi-unit Abutment

Multi-unit Abutment

(Straight)

- Multi-unit Abutment Screw(AANMUS20) included.
- · Use with Multi-unit Driver.
- TCMMUDS20 (short)
- TCMMUDL20 (long)
- Recommend torque: 35Ncm

Cuff Height (mm)	Туре	Ref.C
2.0	Hex	AANMUH5020T
3.0		AANMUH5030T
4.0		AANMUH5040T
5.0		AANMUH5050T
2.0		AANMUN5020T
3.0	Non-Hex	AANMUN5030T
4.0		AANMUN5040T
5.0		AANMUN5050T




Multi-unit Angled Abutment (17°)

- Multi Post Screw(MUMSF/MUMST) included.

Recommend torque: 35Ncm

Cuff Height (mm)	Туре	Ref.C
1.0	Hex Non-Hex	AANMUH50117L
2.0		AANMUH50217L
3.0		AANMUH50317L
4.0		AANMUH50417L
1.0		AANMUN50117L
2.0		AANMUN50217L
3.0		AANMUN50317L
4.0		AANMUN50417L




Multi-unit Angled Abutment (29°)

- Multi Post Screw(MUMSF/MUMST) included.

Recommend torque: 35Ncm

Cuff Height (mm)	Туре	Ref.C
1.0	Hex Non-Hex	AANMUH50129L
2.0		AANMUH50229L
3.0		AANMUH50329L
4.0		AANMUH50429L
1.0		AANMUN50129L
2.0		AANMUN50229L
3.0		AANMUN50329L
4.0		AANMUN50429L



○ Components for Multi-unit Abutment (Continued)

Lab Analog

Profile Diameter	Ref.C	
Ø4.8	RELA300	



Temporary Cylinder

- Cylinder Screw (TASH140) included
- Recommend torque: 15Ncm

Profile Diameter	Туре	Ref.C
Ø4.8	Hex	ETH100T
	Non-Hex	ETN100T



EZ Post Cylinder

- Cylinder Screw (TASH140) included
- Recommend torque : 15Ncm

Profile Diameter	Туре	Ref.C
Ø5.0	Hex	RCA900T
	Non-Hex	RCA800T



Healing Cap

Profile Diameter	Ref.C
Ø5.0	REC600



Impression Coping (Transfer)

Profile Diameter	Ref.C	
Ø4.8	RITE480	



Impression Coping (Pick-Up)

- Guide Pin (RICG150) included

Height (mm)	Туре	Ref.C	
9.4	Hex	RIEH480T	
	Non-Hex	RIEN480T	



Components for Multi-unit Abutment

Gold Cylinder

- Cylinder Screw (TASH140) included
- Useful to make a customized abutment in difficult situations.
- · Precious and non-precious alloys.
- Melting point of gold alloy: 1063°C
- Threaded sleeves for convenient Resin / Wax-up.
- Recommend torque : 15Ncm

Profile Diameter	. Sleeve Color	Ref.C
Ø5.0	Red	REGC200T
	White	REGC100T



CCM Cylinder

- Cylinder Screw (TASH140) included
- Useful to make a customized abutment in difficult situations.
- Can be casted with non-precious alloys (Ni-Cr, Cr-Co alloys).
- Non-precious melting temperature : Depends on Manufacturer
- Threaded sleeves for convenient Resin / Wax-up.
- Melting temperature of CCM : 1300~1400°C
- Recommend torque : 15Ncm

Profile Diameter	Sleeve Color	Ref.C	
Ø4.8	Pink	RCA5013HT	
	Yellow	RCA5013NT	



Plastic Cylinder

- Cylinder Screw (TASH140) included
- Recommend torque : 15Ncm


Profile Diameter	Sleeve Color	Ref.C
Ø5.2	Red	RPEH100T
	White	RPEN100T



II. Abutment Level Prosthesis

4. Flat Abutment & Components

- :The main advantage of this Flat Abutment is the freedom on angulation. Flat Abutment can cover any angulation problems.
- : Only for multiple (Cannot be used for single implant)
- : Only with screw-retained prosthetics.



Components for Flat Abutment (Continued)

Flat Abutment

- Use Hand Driver (1.6 Hex)
- Recommend torque : 25Ncm

Profile Diameter	Cuff Height (mm)	Ref.C
Ø3.5	1	AANFAL3510
	2	AANFAL3520
	3	AANFAL3530
	4	AANFAL3540
	5	AANFAL3550



Flat Cover Screw

• Recommend torque : by hand (5 - 8Ncm)

Profile Diameter	Ref.C
Ø3.5	FCS3510



Flat Healing Abutment

• Recommend torque : by hand (5 - 8Ncm)

Height(mm)	Ref.C
2	FHA402
3	FHA403
4	FHA404



Flat Impression Coping (Transfer)

- Guide Pin (FGPT) included.
- Should be tightened with Impression Driver
- Special impression coping screw which can be used with a 1.2mm hex driver is available on request.


Profile Diameter	Height (mm)	Ref.C	
Ø3.5	9.5	FIT4012T	



Flat Impression Coping (Pick-Up)

- Guide pin (FGPP15) included.

Profile Diameter	Height (mm)	Ref.C	
Ø3.5	12	FIP4012T	



Flat Lab Analog

Profile Diameter	Height (mm)	Ref.C	
Ø3.5	12	FLA3512	



Flat Temporary Cylinder

- Flat Cylinder Screw (FAS) included.
- Recommend torque: 15Ncm

Profile Diameter	Ref.C
Ø4.0	FTC4012T



Flat EZ Post Cylinder

- Flat Cylinder Screw (FAS) included.
- Recommend torque : 25Ncm

Height (mm)	Ref.C
5.5	FEC4005T
7.0	FEC4007T



Flat Gold Cylinder

- Flat Cylinder Screw (FAS) included.
- Useful to make a customized abutment in difficult situations.
- · Precious and non-precious alloys.
- Melting point of gold alloy : 1400 1450°C
- Threaded sleeves for convenient Resin / Wax-up.
- Recommend torque : 25Ncm

Profile Diameter	Ref.C
Ø3.8	FGC4012T



Flat CCM Cylinder

- Flat Cylinder Screw (FAS) included.
- Useful to make a customized abutment in difficult situations.
- Can be casted with non-precious alloys (Ni-Cr, Cr-Co alloys).
- Non-precious melting temperature : Depend on Manufacturer
- Threaded sleeves for convenient Resin / Wax-up.
- Melting temperature of CCM : 1300~1400°C
- Recommend torque : 25Ncm

Profile Diameter	Ref.C
Ø3.8	FCC4012T



Flat Plastic Cylinder

- Flat Cylinder Screw (FAS) included.
- Recommend torque : 25Ncm

Profile Diameter	Ref.C	
Ø4.0	FPC4012T	



NEW PRODUCT

III. Overdenture Prosthesis

1. MegaGen Overdenture System

Meg-Loc

Compatible with products L and K, excellent functionality, & incomparable price!

Combination of Titanium housing and Pekkton (reinforced plastic) creates low water solubility and higher wear resistance and durability than other existing products.

Retention insert offers wide range of retention forces (600gf, 1200gf, 1800gf) to suit each patient, resulting in high level of satisfaction for both patient and dentist. Strong physical properties of Pekkton and insert gap increase elasticity, so that insert does not tear or break unlike conventional nylon products, thereby ensuring strong retention and longer life.





Meg-Ball

Smallest housing, retentive ring with longer life! Even when the implant angle is not parallel, a stable denture can still be produced!

Compatible with other products with Ø2.25 head size, minimized patient inconvenience due to small-size housing, simpler to arrange artificial teeth as space occupied by denture is reduced, and easier to maintain than other systems.

Retentive ring has a high elasticity, abrasion resistance, and durability, thereby doubling the length of life when compared to a silicone O-ring and guaranteeing a longer life than NBR products.

Positioner (0/5/10/15 degrees) maintains parallel housing direction, even with distorted implant placement angle, ensuring denture stability.

Meg-Magnet

Designed to maintain stable and sufficient magnetic force! Completely blocks bursts and corrosion resistant!

Structure is connected with abutment using magnetic force, which is feasible even with insufficient bone volume or poor bone quality

Easy to attach and detach, and minimal inflammation.


Magnet of Ø4.5 & Ø5.0 is compatible with other products, and laser marking on upper part makes it easy to distinguish between up and down.

Sufficient magnetic force (450gf, 650gf) ensures stable retention

Laser sealing blocks any bursting phenomenon.

TiN coating provides corrosion resistance.

Positioner (small & regular) prevents magnet from slipping in the mouth and stops any flow of impression materials under the abutment.





Meg-Rhein

Can compensate for tilted implant placement angle up to 50 $^{\circ}$.


Combined head and housing structure is smallest on the market.

Retentive cap is based on Italian technology and has uniform physical properties. Various retention forces (600gf, 1200gf, 1800gf, 2700gf) classified by color can be selected according to each patient.

Dynamic housing with double structure enables tilting to 25 $^{\circ}$ angle, allowing stable denture even when with distorted implant placement angle.

III. Overdenture Prosthesis

2. Meg-Loc Abutment & Component



▶► Meg-Loc Overdenture System

Advantages

Easy compatibility

Compatible with Product L and Product K (same specifications)

Better abrasion resistance and durability


Combination of Titanium housing and reinforced plastic (Pekkton) provides low water solubility and high resistance, making it superior in abrasion resistance and durability compared to existing products.

Water Sorption Test

Property	Meg-Loc (Pekkton)	Product L	Unit
Water Sorption	8.7	93.5	μg/mm³

Stronger retention and longer life

Strong physical properties of Pekkton and gap in insert increase the elasticity, preventing the insert from being torn or broken unlike existing nylon products, even when angle does not match when attaching & removing denture.



Easy to use

High resistance to plaque and easy cleaning Easy replacement of retention insert



Various Retentive Caps of the Meg-Loc



○ Meg-Loc Overdenture System

Meg-Loc Abutment

- -Angle compensation to one side 20 $^{\circ}$ (both sides 40 $^{\circ})$
- Gently rounded shape
- Compatible with 1.2 Hex Driver
- Recommend torque : 35Ncm

Cuff Height (mm)	Ref.C
0	MLAR00
1.0	MLAR01
2.0	MLAR02
3.0	MLAR03
4.0	MLAR04
5.0	MLAR05
6.0	MLAR06



Meg-Loc Package

- 1 Meg-Loc Abutment
- * Following package items are delivered with San DreMetto Korea packaging.
- 1 Titanium Housing
- 1 Block Out Spacer
- 4 Pekkton Retention Inserts (Gray-250~300gf(for lab), Yellow-600gf, Red-1200gf, Mint-1800gf)

Cuff Height (mm)	Ref.C
0	MLAR00P
1.0	MLAR01P
2.0	MLAR02P
3.0	MLAR03P
4.0	MLAR04P
5.0	MLAR05P
6.0	MLAR06P



Multi Tool

- Retention Insert Insertion & Removal Tool

Ref.C	
MLMT	



III. Overdenture Prosthesis

3. Meg-Ball Abutment & Component



▶► Meg-Ball Overdenture System

Advantages

Easy compatibility



Ø2.25 head size for easy compatibility with other products

Smallest Housing



Small housing minimizes patient inconvenience, facilitates arrangement of artificial teeth by reducing space occupied by denture, and is easier to maintain than other systems.

Double length of life



High elasticity, abrasion resistance, and durability doubles the length of life when compared with silicone O-ring and guarantees longer life than NBR products.

Stable denture even when implant placement angle is distorted

Positioner (0/5/10/15 degrees) maintains parallel housing direction even when angle of implant placement is distorted, ensuring denture stability









Tilting Angle



○ Meg-Ball Overdenture System

Meg-Ball Abutment

- Angle compensation to one side 15 $^{\circ}$ (both sides 30 $^{\circ})$
- Ø2.25 Ball shape
- Recommend torque : 35Ncm

Cuff Height (mm)	Ref.C
0	MBAR00
1.0	MBAR10
2.0	MBAR20
3.0	MBAR30
4.0	MBAR40
5.0	MBAR50
6.0	MBAR60



Meg-Ball Package

- Composed of Meg-Ball Abutment/ Metal Housing Set/ Housing Positioner (0°,5°,10°,15°)

Cuff Height (mm)	Ref.C
0	MBAR00P
1.0	MBAR10P
2.0	MBAR20P
3.0	MBAR30P
4.0	MBAR40P
5.0	MBAR50P
6.0	MBAR60P



Meg-Ball Metal Housing Set

- 1 Metal Housing
- 1 Retentive Ring

Ref.C	
MBHR	




Retentive Ring Set

Ref.C
MBR5
MBR10



III. Overdenture Prosthesis

4. Meg-Magnet Abutment & Component



Meg-Magnet Overdenture System

Advantages

Easy to apply for elderly patients or disabled patients

Applicable with insufficient bone volume and poor bone quality Easy to attach and detach Unlikely to cause inflammation

Designed for maximum magnetic efficiency and durability

Sufficient magnetic force (450gf, 650gf) to ensure stable retention Laser sealing blocks any bursting phenomenon

Outstanding retention

- Blocks bursting
- Corrosion resistant
- Abrasion resistant

TiN coating provides corrosion resistance Over 0.1mm thickness at contact with attachment to ensure wear resistance



Easy to distinguish between up and down via laser marking on upper section

Magnet of Ø4.5 & Ø5.0 is compatible with other products Laser marking on upper part makes it easy to distinguish between up and down





No slippage of magnet

Positioner (small & regular) prevents magnet from slipping in mouth and stops any flow of impression materials under the abutment





(450af)

Ø5.0

Component of the Meg-Magnet




○ Meg-Magnet Overdenture System

Meg-Magnet Abutment

- Use to 1.2 Hex Driver
- Recommend torque : 35Ncm

Profile Diameter	Cuff Height (mm)	Ref.C
	0	MMAR400
	1.0	MMAR410
04.5	2.0	MMAR420
Ø4.5	3.0	MMAR430
	4.0	MMAR440
	5.0	MMAR450
	0	MMAR500
	1.0	MMAR510
Ø5.0	2.0	MMAR520
Ø5.0	3.0	MMAR530
	4.0	MMAR540
	5.0	MMAR550



Meg-Magnet Package

- 1 Meg-Magnet Abutment
- 1 Magnet (Ø4.5-450gf, Ø5.0-650gf)
- 1 Magnetic Positioner

Profile Diameter	Cuff Height (mm)	Ref.C
	0	MMAR400P
	1.0	MMAR410P
Ø4.5	2.0	MMAR420P
<i>1</i> 04.5	3.0	MMAR430P
	4.0	MMAR440P
	5.0	MMAR450P
	0	MMAR500P
	1.0	MMAR510P
Ø5.0	2.0	MMAR520P
25.0	3.0	MMAR530P
	4.0	MMAR540P
	5.0	MMAR550P




Meg-Magnet Attachment Set

Size	Ref.C
Small	MA402
Regular	MA502



III. Overdenture Prosthesis

5. Meg-Rhein Abutment & Components



Meg-Rhein Overdenture System

Meg-Rhein Overdenture System

(Dynamic)

- 1 Meg-Rhein Abutment
- 1 Plastic Impression Coping
- 1 Stainless Steel Housing (Dynamic) & Black-Lab
- 1 Protective Disk
- 4 Retentive Caps (Yellow-0.6kgf, Pink-1.2kgf, White-1.8kgf, Violet-2.7kgf)
- · Perfect compatibility with the Rhein83 from Italy.
- · Recommend torque: 15Ncm.





Meg-Rhein Abutment with Plastic Impression Coping

Meg-Rhein Overdenture System


- 1 Meg-Rhein Abutment
- 1 Plastic Impression Coping
- 1 Stainless Steel Housing
- 1 Protective Disk
- 5 Retentive Caps (Black-Lab, Yellow-0.6kgf, Pink-1.2kgf, White-1.8kgf, Violet-2.7kgf)
- · Perfect compatibility with the Rhein83 from Italy.
- Recommend torque: 15Ncm.

Cuff Height (mm)	Ref.C
0	ADR00P
1.0	ADR01P
2.0	ADR02P
3.0	ADR03P
4.0	ADR04P
5.0	ADR05P
6.0	ADR06P



with Plastic Impression Coping

▶▶ Overdenture System



Components for Meg-Rhein Abutment (Continued)

Stainless Steel Housing

• 5ea/pack

Ref.C	
MHP	



Stainless Steel Housing

(Dynamic)

• 5ea/pack

Ref.C	
THP	



Retentive Caps (White)

- White cap(1.8kg) For refill (5ea/pack).
- Can be used for more retentive force following pink cap(1.2kgf).

Ref.C)
RCWP	



Retentive Caps (Violet)

- Violet cap(2.7kg) For refill (5ea/pack).
- Can be used for more retentive force following white cap(1.8kgf).

Ref.C
RCVP



Retentive Caps (Pink)

• Pink cap(1.2kgf) - For refill (5ea/pack).


Ref.C	
RCPP	



Components for Meg-Rhein Abutment

Retentive Caps (Yellow)

• Yellow cap(0.6kgf) - For refill (5ea/pack).





Retentive Caps (Black)

• For laboratory

Ref.C
RCBP



Stainless Impression Coping (Pick-Up)


- · 2ea/pack.
- Italy Rhein 83 products.
- For accurate (pick-up type) impression.
- Metal with groove design to prevent from swaying.

Ref.C	
044CAIN	





Lab Analog






Retentive Cap Removal Tool

Retentive Cap removal tool.

Ref.C 091EC



Retentive Cap Insertion Tool

Retentive Cap insertion tool.

Ref.C 085IAC



AnyRidge Kit


I. AnyRidge Abutment Selection Guide Kit

Ref.C KANASG3000

For the best selection of abutments.

- Colors indicate different cuff heights (Yellow: 2mm, White: 3mm, Pink: 4mm, Red: 5mm).
- · Store 2 pieces in each container.
- Autoclavable to sterilize.


















Angle type(15°)
(Angled Abutment select)

Angle type(25°)
(Angled Abutment select)

(EZ Post & Solid Abutment select)

II. AnyRidge Surgical Kit : Standard Type


Ref.C KARIN3003



II. AnyRidge Surgical Kit : Full Type



Easier and safer to drill for the depth as you need with the stopper drills.



Surgical Kit Components

Lance Drill

• Useful to make an indentation on cortical bone to confirm the exact drilling location.

Diameter	Туре	Ref.C	
Ø2.0	Long	MGD100L	



Marking Drill


- Each drill has calibrations from 7.0 to 18.0mm. (TANSDF4815, TANSDF5415, TANSDF5915 have calibrations up to 15.0mm)
- Easy to recognize by dual marking systems. (Groove and laser marking)

Diameter	Length (mm)	Ref.C
Ø2.0		TANTDF2018
Ø2.5		SD2518S
Ø2.8	10	SD2818S
Ø3.3	18	TANSDF3318
Ø3.8		TANSDF3818
Ø4.3		TANSDF4318
Ø4.8	15	TANSDF4815
Ø5.4		TANSDF5415
Ø5.9		TANSDF5915



Stopper Drill

Diameter	Length (mm)	Ref.C
	7	TANTDF2007
0 0.0	8.5	TANTDF2008
Ø2.0	10	TANTDF2010
	11.5	TANTDF2011
	7	SD2807M
0 0.0	8.5	SD2808M
Ø2.8	10	SD2810M
	11.5	SD2811M
	7	TANSDF3307
20.0	8.5	TANSDF3308
Ø3.3	10	TANSDF3310
	11.5	TANSDF3311
	7	TANSDF3807
0 0.0	8.5	TANSDF3808
Ø3.8	10	TANSDF3810
	11.5	TANSDF3811
	7	TANSDF4307
Q4.0	8.5	TANSDF4308
Ø4.3	10	TANSDF4310
	11.5	TANSDF4311
	7	TANSDF4807
Ø4.0	8.5	TANSDF4808
Ø4.8	10	TANSDF4810
	11.5	TANSDF4811



Surgical Kit Components (Continued)

Point Trephine Bur


Diameter	Ref.C
Ø5.0 (ln.Ø4.0)	SPTB4050



Trephine Bur

- Minimizes the drilling steps needed, especially for wider fixtures.
- Helpful for collecting autogenous bone.
- Useful for removing failed and fractured fixtures.
- Depth markings are 7, 8.5, 10, 11.5, 13mm, same depths as fixtures. (No Y dimension so markings are actual length).
- Markings on the drill shaft represent the inside / outside diameter of Trephine Burs.

Diameter	Туре	Ref.C
Ø3.5 (in Ø2.5)		TANTBL2535
Ø5.0 (in Ø4.0)	Short	TANTBL4050
Ø6.0 (in Ø5.0)	(32mm)	*TANTBL5060
Ø7.0 (in Ø6.0)		*TANTBL6070
Ø3.5 (in Ø2.5)		*TANTBE2535
Ø5.0 (in Ø4.0)	Long	*TANTBE4050
Ø6.0 (in Ø5.0)	(38mm)	*TANTBE5060
Ø7.0 (in Ø6.0)		*TANTBE6070




- (*) Separate sales item.
- 3.5, 5.0 Trephine Bur are included in Surgical kit.

Cortical Bone Drill

- Removes cortical bone and enlarges osteotomy socket especially at hard bone.
- Similar function with countersink drill of other systems.
- Each drill has two steps of diameter for convenience.

Diameter	Ref.C
Ø3.5	TANCDL3500
Ø4.0~ Ø5.5	TANCDL4055
Ø6.0~ Ø8.0	TANCDL6080





Handpiece Connector

- Delivers torque for the placement of a fixture with a handpiece.
- · Easy and secure pick-up and delivery.
- · Used to place an implant without a mount.
- Marks on the shaft can indicate the position of fixture platform, especially in flapless surgery.

Length (mm)	Туре	Ref.C
5	*Ultra short	TANHCU
10	Short	TANHCS
15	Long	TANHCL
10	Short (MiNi)	HCS17
15	Long (MiNi)	HCL17

(*) Separate sales item.





Ratchet Connector

- Delivers torque for the placement or removal of a fixture with a Ratchet Wrench.
- Secure a Ratchet Extension or Torque Wrench to a fixture before exerting force.
- Too much torque force can result a damage to the hex of a fixture.
- Marks on the shaft can indicate the position of fixture platform, especially for flapless surgery.

Length (mm)	Туре	Ref.C
6	*Ultra short	TANREU
10	Short	TANRES
15	Long	TANREL
15	Short(MiNi)	RCS17
20	Long (MiNi)	RCL17

(*) Separate sales item.





Hand Driver (1.2 Hex)

- Used for all Cover Screws, all Abutment Screws and all Healing Abutments.
- · Available in 4 lengths for convenience.
- Hand Driver can be directly inserted into the Torque Wrench without using an adapter.
- Hex tip can withstand 35-45Ncm of torque without distortion.

Length(mm)	Туре	Ref.C
5	*Ultra-short	TCMHDU1200
10	Short	TCMHDS1200
15	Long	TCMHDL1200
20	*Extra-long	TCMHDE1200

(*) Separate sales item.



Surgical Kit Components

Abutment Removal Driver

- Used to remove final abutment; use after removing Abutment Screw.
- · Insert straight into the abutment and rotate clockwise.
- Long Abutment Removal Driver is for disconnecting an abutment with a cemented crown.


Length (mm)	Ref.C
17.5	TANMRD18
25.0	*TANMRD25

(*) Separate sales item.



Drill Extension

- Extends drills & other handpiece tools.
- No more than 35Ncm torque: Can be distorted when too much force is applied.






Direction Indicator

- Confirms drilling direction and location during drilling.
- · Checks drilling position.

Length (mm)	Ref.C
Ø2.0 / Ø2.8	MDI100
Ø3.2 / Ø4.7	MDI3348
	Ø2.0 / Ø2.8



Path Finder

- After placing a fixture, a Path Finder can be connected to guide parallel for the next implant.
- Gingival depth can be measured with the grooves especially for flapless surgeries.

Length (mm)	Ref.C
10	TANPFF3580




Torque Wrench & Adapter

 Torque Wrench has torque options from 15Ncm to 45Ncm and is used for the placement of an implant and final tightening of the Abutment Screw.

Туре	Ref.C
Torque Wrench	MTW300AT
*Torque Wrench Adapter(Handpiece)	TTAI100
Torque Wrench Adapter(Ratchet)	TTAR100


(*) Separate sales item.



III. AnyRidge Prosthetic Kit

Ref.C KANPK3000

A kit with all kinds of driver that are needed for prosthetics.



Prosthetic Kit Components

Solid Driver

- For the delivery of Solid Abutments.
- Color coded for different profile diameters.
 (Ø4-magenta, Ø5-blue, Ø6-yellow, Ø7-green)
- Two different heights. (8.5 / 13.5mm)
- Directly connectable to Torque Wrench.

Solid Abutment Profile Diameter	Length(mm)	Ref.C	
Ø4	8.5	TANSDS400	
<i>Ø</i> 4	13.5	*TANSDL400	
Q.F.	8.5	TANSDS500	
Ø5	13.5	*TANSDL500	
90	8.5	TANSDS600	
Ø6	13.5	*TANSDL600	
Ø7	8.5	TANSDS700	
	13.5	*TANSDL700	



(*) Separate sales item.

Octa Driver

- For seating of the Octa Abutment into the fixture.
- Can also be connected to Torque Wrench.

Length (mm)	Ref.C
7	MOD300S
13	MOD300L




Ball Driver

- For seating of the Ball Abutment into the fixture.
- Can connect to a Handpiece, Ratchet or Torque Wrench.
- · Available in long and short.

Туре	Ref.c
*Handpiece Connector(Short)	TBH250S
*Handpiece Connector(Long)	TBH250L
*Ratchet Connector(Short)	TBR250S
*Ratchet Connector(Long)	TBR250L
Toque Driver(Short)	TBT250S
*Toque Driver(Long)	TBT250L


(*) Separate sales item.



Impression Coping Driver (Transfer)

- For transfer type of Impression Coping.
- Works with friction only.
- Small but powerful grip.

Туре	Ref.C
For Two piece impression Coping	TCMID
For One piece impression Coping	TCMIDE




- 1. Connect Impression Coping and Impression Driver together
- 2. Adjust Connection with a Fixture by turning a Holder clockwise.
- 3. Push the Holder and put the Impression Coping into the Fixture.
- 4. Turn the Driver clockwise to ensure connection of the Impression Coping and Fixture.

IV. AnyRidge Bone Profiler Kit

Ref.C KARBP3000

Removes the overhanged bone around a fixture to allow adequate seating of a Healing Abutment or a Prosthetic Abutment.

- Place a Guide Pin into a fixture and choose a Bone Profiler which fits with the situation.
- Four different sizes of bone profiler and four guide pins are included in the kit.



Bone Profiler

- Guide Pin(TANPGF3305) included.
- Each bone profiler can be purchased separately for refill.
- Each pakage includes a bone profiler and a guide pin.

Profile Diameter	Length (mm)	Ref.C	
Ø4	13	TANBPL40G	
Ø5		TANBPL50G	
Ø6	8	TANBPS60G	
Ø7		TANBPS70G	



V. Optional components

- not included in the surgical kit
- can be purchased separately and placed into the 'option' spaces provided in the surgical kit

Right Angle Driver Tip

- Used for all Cover Screws, all abutment screws and all Healing Abutments.
- Hex tip can withstand 35-45Ncm of torque without distorting.

Length(mm)	Туре		Ref.C
4	Ultra-short		MDR120SS
10	Short	Hex 1.2	MDR120S
15	Long		MDR120L
20	Extra Long		MDR120EL



Lindermann Drill

- · Cross cut on the drill.
- · Can correct the path during drilling.

Diameter(mm)	Ref.C
2	TEEL200M



Insert Driver

- Used for all Cover Screws, all abutment screws and all Healing Abutments.
- Hex tip can withstand 35-45Ncm of torque without distorting.

Length(mi	m) T	ype	Ref.C	
10	Short	Llau 4.0	MID120S	
15	Long	Hex 1.2	MID120L	



Hand Tap

- Useful when the internal screw of a fixture is damaged.
- Retapping damaged threads.
- Need to be patient and force-controlled.

Туре	Ref.C
M1.8	THT180L



Multi-unit Driver (2.0 Hex) (For Multi-unit Abutment)

• For the seating & tightening of multi-unit Abutment (Straight type)

Length(mm)	Туре	Ref.C	
10	Short	TCMMUDS20	
15	Long	TCMMUDL20	



Flattening Drill

- In the case of irregular bone, stopper drill can be drilled in precise depth by flattening the bone.
- Flattening Lance and Housing are connected together. Two types of Housing diameters (Ø5.0 & Ø6.0) are composed in accordance with the size of final drill diameter.
- Ø5.0 = Stopper Drill Ø2.0~ Ø4.3
- Ø6.0 = Stopper Drill Ø4.8~ Ø5.4
- Formation of boundary through housing will guide the next drilling location of fixture.

Diameter	Length (mm)	Ref.C
Ø5.0 / Ø2.0	0.5	FD5020
Ø6.0 / Ø2.0	3.5	FD6020



Manual Inserter

- Specially designed for manual placement of AnyRidge fixture.
- Especially useful at immediate implant placement on maxillary anterior.
- The tip has same structure with the hand-piece connector.

Ref.C	
TANMI	



Reamer Drill & Center Pin

- Removes inner lip of the cast after casting Burn-out Cylinders of Solid Abutment.
- Center Pin have 4 different diameters according to the profile diameter of Solid Abutments.


Diameter	Туре	Ref.C
Ø10.0	Reamer Drill	TANRD
Ø4.0	Center Pin	TANRDJ40
Ø5.0		TANRDJ50
Ø6.0		TANRDJ60
Ø7.0		TANRDJ70



Trephine Bur Stopper

- Controls the depth of trephination with a Stopper placed into the Trephine.
- Especially useful in cases with limited availabe bone from important anatomy.


Length (mm)	Ref.C
7.0	TANTSF2307
8.5	TANTSF2308
10.0	TANTSF2310
11.5	TANTSF2311



Bottom Drill


- It removes remaining bone in osteotomy socket after trephine drilling.
- It imprints the sizes of fixtures, for example 7, 8.5, 10, 11.5 and 13mm, by laser marker.

Diameter	Туре	Ref.C
Ø3.3		TCMBDS33
Ø3.8		TCMBDS38
Ø4.8	Short (32mm)	TCMBDS48
Ø5.8	(OZITIITI)	TCMBDS58
Ø6.8		TCMBDS68
Ø3.3		TCMBDL33
Ø3.8		TCMBDL38
Ø4.8	Long (38mm)	TCMBDL48
Ø5.8		TCMBDL58
Ø6.8		TCMBDL68



Ratchet Wrench

- Used to exert more force than handpiece.
- No bearing system: No breakage and corrosion problems.
- Attaches to Ratchet Extension.
- Arrow laser marking indicates direction of force.







R2GATE Full Surgical KIT


I. R2GATE Full Surgical Kit for AnyRidge System

- If you only use a specific system, corresponding system's full kit can be provided.
- R2GATE full surgical kit is composed with all of drills and components that are needed for the Digital Guided Surgery which uses R2GATE Guide™ after R2GATE™ diagnosis. It helps to actualize minimally invasive surgery and makes exact clinical result as the diagnosis.









Components for R2GATE Full Surgical Kit (Continued)

- If you only use a specific system, corresponding system's full kit can be provided.
- R2GATE full surgical kit is composed with all of drills and components that are needed for the Digital Guided Surgery which uses R2GATE GuideTM after R2GATETM diagnosis. It helps to actualize minimally invasive surgery and makes exact clinical result as the diagnosis.

Initial Drill

- Use the initial drill in order to mark the drilling position on the bone. Start drilling slowly, when drill guide part is fully contacted with drilling core of R2GATE Guide™.
- Recommended drilling speed range is 300 ~ 800 RPM with copious irrigation.


Diameter	Guide Diameter	Length(mm)	Ref.C	
Ø2.6	Ø5.0	1.0	R2ID2601	



Second Drill

- This unique step-drill(from Ø2.0 to Ø4.6) is used to flare out the upper cortical bone of the osseotomy.
- It helps not only the rest drilling procedure but abut- ment connection. In case of hard bone, if the 2nd drilling will be disturbed by thick cortical bone. Stop the drilling and try it after final drilling procedure.

Diameter	Guide Diameter	Length(mm)	Ref.C
Ø2.5	Ø5.0	5.0	R2SD2505



Stopper Drill

- Universal drills consist of Ø2.0, Ø.2.5, Ø2.8 diameter to enlarge the osteotomy gradually.
- The length of drill are designed as 7.0, 8.5, 10, 11.5,13mm for most common length of implant system.
- Recommended drilling speed range is 500 \sim 800 RPM with copious irrigation.

Diameter	Guide Diameter	Length(mm)	Ref.C
		6.5	AGSD2007
		8.0	AGSD2008
Ø2.0		9.5	AGSD2010
		11.0	AGSD2011
		12.5	AGSD2013
	Ø5.0	6.5	AGSD2507
		8.0	AGSD2508
Ø2.5		9.5	AGSD2510
		11.0	AGSD2511
		12.5	AGSD2513
		6.5	AGSD2807
		8.0	AGSD2808
Ø2.8		9.5	AGSD2810
		11.0	AGSD2811
		12.5	AGSD2813



Bone Profiler


• Recommended drilling speed is 300 ~ 800 RPM.

Diameter	Guide Diameter	Ref.C
Ø4.0	Ø5.0	AGBP40
Ø5.0		AGBP50
Ø6.0	Ø6.5	AGBP60



Stopper Drill

• Recommended drilling speed is 300 ~ 800 RPM.




	Diameter	Guide Diameter	Length(mm)	Ref.C
			6.5	ARSD3307
			8.0	ARSD3308
	Ø3.3		9.5	ARSD3310
			11.0	ARSD3311
		Ø5.0	12.5	ARSD3313
	Ø3.8		6.5	ARSD3807
			8.0	ARSD3808
			9.5	ARSD3810
			11.0	ARSD3811
			12.5	ARSD3813
			6.5	ARSD4307
			8.0	ARSD4308
	Ø4.3		9.5	ARSD4310
			11.0	ARSD4311
			12.5	ARSD4313

Diam	eter	Guide Diameter	Length(mm)	Ref.C
			6.5	ARSD4807
			8.0	ARSD4808
Ø4.8	8		9.5	ARSD4810
			11.0	ARSD4811
		Ø6.5	12.5	ARSD4813
			6.5	ARSD5407
			8.0	ARSD5408
Ø5.	4		9.5	ARSD5410
			11.0	ARSD5411
			12.5	ARSD5413
			6.5	ARSD5908
			8.0	ARSD5907
Ø5.	9		9.5	ARSD5910
			11.0	ARSD5911
			12.5	ARSD5913

Cortical Bone Drill

• Recommended drilling speed : 300 ~ 800 RPM

Diameter	Guide Diameter	Length(mm)	Ref.C
Ø3.4		5.0	R2CD3405
Ø3.8	ØF 0		R2CD3805
Ø4.3	Ø5.0 Ø6.5		R2CD4305
Ø4.8			R2CD4805
Ø5.3			R2CD5305
Ø5.8			R2CD5805
Ø6.3			R2CD6305



Components for R2GATE Full Surgical Kit

Implant Carrier

- The purpose of tab drills in the universal kit system is insertion test. some of implant are required this procedure before final fixture insertion. choose the one-step under size of tab to protect from enlargement of osteotomy.
- Recommended insertion torque and speed is 45 \sim 50Ncm, under 40 RPM.

Connection	Guide Diameter	Туре	Ref.C
	Ø5.0	Ratchet	ICRH2324
2.3 Hex	Ø6.5		ICWH2324
2.3 HeX	Ø5.0	Handpiece	ICRH2324H
	Ø6.5		ICWH2324H



Carrier-Handpiece Adapter

• Useful to use the handpiece for the implant placement following initial delivery of a fixture with a fixture carrier.

Diameter	Ref.C
5.0	AGHA



Carrier Extension


• To extend the length of implant carrier.

Diameter	Ref.C
4.0	MRE400S



Drill Extension

- No more than 35Ncm torque : May distorted when excessive force is applied.
- Extends drills & other handpiece instruments.





Hand Driver (1.2 Hex)

- Used for all Cover Screws, Abutment Screws, and Healing Abutments.
- Available in 4 lengths for added convenience.
- Hand Driver can be directly inserted into the Torque Wrench without using an adaptor.
- Hex tip can with stand 35-45Ncm of torque without distorting.

Length(mm)	Туре	Ref.C
5.0	*Ultra-short	TCMHDU1200
10	Short	TCMHDS1200
15	Long	TCMHDL1200
20	*Extra-long	TCMHDE1200





Ratchet Wrench

- ${\boldsymbol{\cdot}}$ Used to exert more force than the Handpiece.
- No bearing system : No breakage and no corrosion problems.
- Arrow laser marking indicates direction of force.






R2GATE Universal Kit

Maximize the cost-effectiveness & efficiency.


Ref.C KAGUN3000

When you want to do R2GATE surgery with R2GATE Guide™, Please inform us your favorite implant system

Make your own R2GATE Surgical Kit with your favorite implant system. R2GATE Universal kit consists of basic drilling set which can be used for any implant system. You can add system options as "Implant Carrier", "Cortical Bone Drill", "Tap Drill" to your favorite implant system. The specification of final drills will be decided with treatment planning and delivered to you with R2GATE Guide™ will be from the R2GATE Design Center.




R2GATE Universal Kit



Drills & Components for R2GATE Universal Kit

Basic drilling set for any implant system. It consists of initial drill, 2nd drill, universal drills and essential tools.



Initial Drill

- Use the initial drill in order to mark the drilling position on the bone. Start drilling slowly, when drill guide part is fully contacted with drilling core of R2GATE GuideTM.
- Recommended drilling speed range is 300 ~ 800 RPM with copious irrigation.

Diameter	Guide Diameter	Length(mm)	Ref.C	
Ø2.6	Ø5.0	1.0	R2ID2601	

16

12.5

2.25

Second Drill

- This unique step-drill(from $\varnothing 2.0$ to $\varnothing 4.6$) is used to flare out the upper cortical bone of the osseotomy.
- It helps not only the rest drilling procedure but abutment connection. In case of hard bone, if the 2nd drilling will be disturbed by thick cortical bone. Stop the drilling and try it after final drilling procedure.

Diameter Guide Diameter		Length(mm)	Ref.C	
Ø2.5	Ø5.0	5.0	R2SD2505	



Stopper Drill

- Universal drills consist of @2.0, @.2.5, @2.8 diameter to enlarge the osteotomy gradually.
- The length of drill are designed as 7.0, 8.5, 10, 11.5,13mm for most common length of implant system.
- Recommended drilling speed range is 500 ~ 800 RPM with copious irrigation.

Diameter	Guide Diameter	Length(mm)	Ref.C
		6.5	R2SD2007
		8.0	R2SD2008
Ø2.0		9.5	R2SD2010
		11.0	R2SD2011
		12.5	R2SD2013
	Ø2.5 Ø5.0	6.5	R2SD2507
		8.0	R2SD2508
Ø2.5		9.5	R2SD2510
		11.0	R2SD2511
		12.5	R2SD2513
		6.5	R2SD2807
		8.0	R2SD2808
Ø2.8		9.5	R2SD2810
		11.0	R2SD2811
		12.5	R2SD2813



Carrier-Handpiece Adapter

 Useful to use the handpiece for the implant placement following initial delivery of a fixture with a fixture carrier ratchet type.

Square	Ref.C
4.0	AGHA



Carrier Extension

To extend the length of implant carrier.

Square	Ref.C
4.0	MRE400S



Torque Wrench & Adapter

 Torque Wrench has torque options from 15Ncm to 45Ncm and is used for the placement of an implant and final tightening of the Abutment Screw.

Туре	Ref.C
Torque Wrench	TW70
Torque Wrench Adapter(Ratchet)	TTAR100






Cortical Bone Drill[AR]

• Recommended drilling speed: 300 ~ 800 RPM


Diameter	Guide Diameter	Length(mm)	Ref.C
Ø3.4	Ø5.0		R2CD3405
Ø3.8			R2CD3805
Ø4.3			R2CD4305
Ø4.8		5.0	R2CD4805
Ø5.3			R2CD5305
Ø5.8	Ø6.5		R2CD5805
Ø6.3			R2CD6305




Optional Instrument

Narrow Crest Drill

- It is used when fixture will be slantly implanted or to flat the sloped bone surface of narrow ridge to prevent any slips during drilling.
- · Design as 2-piece: drill body and housing
- Can be disassembled. Easy to clean and remove bone chips
- Can harvest autogenous bone if it is used after soft tissue



Diameter Guide Diamet		Length(mm)	Ref.C	
Ø4.0	Ø5.0	15.5(12.5/2)	NCD402	





Set the site by drilling counter-clockwisely with low speed (\leq 100rpm)



Start drilling clockwisely (400~600rpm)



Bone is now flat. Perform drilling with proper drilling sequence.






Disassemble body and housing after drilling to remove bone chip. Clean and sterilize after every usage.

○ Final Drill Option [Disposable or Permanent]

Stopper Drill[Straight] For all implant system


- · Common use
- Step back type drillling
- Provided from local R2GATE Design Center to users. The size of disposable drills are decided depend size on treatment planning regarding to fixture size and bone density of patient.
- Recommended drilling speed is 300 ~ 800
- · Final drill.
- $\boldsymbol{\cdot}\,$ The base is disposable and can be made for permanent under your order


Diameter	Diameter	Length(mm)	Ref.C	Ref.C
		7.0	R2PS3407	R2DS3407
		8.0	R2PS3408	R2DS3408
		9.0	R2PS3409	R2DS3409
Ø3.4		10.0	R2PS3410	R2DS3410
		11.0	R2PS3411	R2DS3411
		12.0	R2PS3412	R2DS3412
		13.0	R2PS3413	R2DS3413
		7.0	R2PS3807	R2DS3807
		8.0	R2PS3808	R2DS3808
		9.0	R2PS3809	R2DS3809
Ø3.8	Ø5.0	10.0	R2PS3810	R2DS3810
		11.0	R2PS3811	R2DS3811
		12.0	R2PS3812	R2DS3812
		13.0	R2PS3813	R2DS3813
		7.0	R2PS4307	R2DS4307
		8.0	R2PS4308	R2DS4308
		9.0	R2PS4309	R2DS4309
Ø4.3		10.0	R2PS4310	R2DS4310
		11.0	R2PS4311	R2DS4311
		12.0	R2PS4312	R2DS4312
		13.0	R2PS4313	R2DS4313
		7.0	R2PS4807	R2DS4807
		8.0	R2PS4808	R2DS4808
		9.0	R2PS4809	R2DS4809
Ø4.8		10.0	R2PS4810	R2DS4810
		11.0	R2PS4811	R2DS4811
		12.0	R2PS4812	R2DS4812
		13.0	R2PS4813	R2DS4813
		7.0	R2PS5307	R2DS5307
		8.0	R2PS5308	R2DS5308
		9.0	R2PS5309	R2DS5309
Ø5.3	Ø6.5	10.0	R2PS5310	R2DS5310
		11.0	R2PS5311	R2DS5311
		12.0	R2PS5312	R2DS5312
		13.0	R2PS5313	R2DS5313
		7.0	R2PS5807	R2DS5807
		8.0	R2PS5808	R2DS5808
		9.0	R2PS5809	R2DS5809
Ø5.8		10.0	R2PS5810	R2DS5810
		11.0	R2PS5811	R2DS5811
		12.0	R2PS5812	R2DS5812
		13/0	R2PS5813	R2DS5813



Drill position on the kit

- Every disposable drills have the numbering system to clarify it's own position on the universal kit.
 Check the drill size and position number, then
- install it to the right position.





Sterilized package

- All disposable drills are packaged at clean room and sterilized by "Gamma-ray".
 Check the "Sterilized" seal on the package and
- open it at the operation site before surgery.



Digital Material

I. ZrGEN®

ZrGEN° is the brand name of MegaGen Titanium Base. ZrGEN provides an aesthetic outcome and simplified dental implant prosthesis. A ZrGEN° crown and monolithic crown connected to a ZrGEN° Abutment provide strong and precise connection with the implant fixture.

Variety of ZrGEN®



ZrGEN Coping



PMMA Provisional Crown



ZrGEN Monolithic



ZrGEN Crown



ZrGEN Bridge



ZrGEN Coping for PFZ

ZrGEN° Sub Structure



ZrGEN



Zirconia customized body



Zirconia Final Crown

ZrGEN°

The strength of ZrGEN° frees you from the chipping of conventional PFM prosthesis. Monolithic zirconia crowns have no metal substructure, ensuring more aesthetic results. ZrGEN° crown and bridge are a superior substitutes for all conventional dental materials.



Tooth shade cuff area



Minimized Ti-connection

Clinical Application













II. TIGEN®

TiGEN° is the brand name of MegaGen Titanium customized abutment. It promises outstanding durability and simplified dental implant prosthesis. Ready-made connection part provides a strong and precise connection with the implant fixture.








○ ZrGEN° Prosthesis



ZrGEN® Abutment

ZrGEN* Abutment provides a strong and precise connection with the implant fixture. With Zirconia* coping, crown margins can be placed supragingivally since zirconia material matches with the color of natural teeth. Residual cement problems are no longer an issue.



Scan Abutment Option

Scan Abutment

- Abutment Screw included. . AnyRidge (AANMSF)
- . AnyRidge Octa 1 (AROAS16B/ AROAS16)
- . AnyOne Internal (AS20)
- . AnyOne Exeternal (SCS160/ RCS200)
- . AnyOne OneStage (EXIMS100)
- MiNi (MIAS14)
- ST (OSGSAS3110/ OSGSAS3210)
- . Octa Level (IRCS200)
- . Multi-unit Abutment (MUAS)
- · For Chairside/ Labside
- · Included spare Abutment Screw
- Surpporting Dental CAD
- 3 Shape
- Exocad
- Dental Wings

System	Profile Diameter	Length (mm)	Туре	Ref.C
AnyRidge	Ø4.0	9	-	AANISR4009T
		13	-	AANISR4013T
AnyRidge	Ø4.0	13	NC	AROSANT
Octa 1			RC	AROSART
AnyOne Internal	Ø4.0	9	-	AAOISR4009T
		13	-	AAOISR4013T
AnyOne External		9	Small	AEXESS4009T
		13		AEXESS4013T
		9	Regular	AEXESR4009T
		13		AEXESR4013T
AnyOne OneStage	Ø4.0	13	Cuff 1.8	AEXISR4010T
MiNi	3.5	9	-	MISS3509T
		13	-	MISS3513T
ST	Ø4.0	9	Small	OSGSSC3110T
		13	Small	OSGSSC3111T
		9	Dogudor	OSGSSC3210T
		13	Regular	OSGSSC3211T
Octa Level	Ø4.0	11	-	AOCESC4011T
MUA Level (N-Type)	Ø4.0	13	-	AMUASR4013T



Scan Abutmet (C-type)

- Abutment Screw included. . AnyRidge (AANMSF) . AnyOne (AS20)
- . AnyRidge Octa 1(AROAS16B/ AROAS16)
- $\bullet \ \, \text{Scan Post for Sirona Cerec users} \to \text{CEREC} \\ \bullet \ \, \text{In in Lab CAD Software, compatible with}$
- Xive Library

Syste	m	Profile Diameter	Cuff Height	Post Size	Ref.C
			0.5		ARICSS3405T
AnyRidge		Ø3.9	1	Small	ARICSS3410T
			2		ARICSS3420T
		Ø4.3	0.5		ARICSS3805T
			1		ARICSS3810T
			2		ARICSS3820T
		Ø5.5	0.5	Large	ARICSL4505T
			1		ARICSL4510T
			2		ARICSL4520T
	NC :	Ø3.9	0.5	Small	AROCSS3405NT
			1		AROCSS3410NT
			2		AROCSS3420NT
		Ø4.3	0.5	Small	AROCSS3805NT
			1		AROCSS3810NT
			2		AROCSS3820NT
A Di-l	RC	Ø3.9	0.5	Small	AROCSS3405RT
AnyRidge Octa 1			1		AROCSS3410RT
Ocia i			2		AROCSS3420RT
		Ø4.3	0.5		AROCSS3805RT
			1		AROCSS3810RT
			2		AROCSS3820RT
		Ø5.5	0.5	Large	AROCSL4505RT
			1		AROCSL4510RT
			2		AROCSL4520RT
			0.5		AOICSS3405T
		Ø3.9	1		AOICSS3410T
AnyOne			2	Small	AOICSS3420T
		Ø4.3	0.5		AOICSS3805T
			1		AOICSS3810T
			2		AOICSS3820T
		Ø5.5	0.5	Large	AOICSL4505T
			1		AOICSL4510T
			2		AOICSL4520T





Scan Healing **Abutment** & Scan Post

- Abutment Screw included. · AnyRidge (ARIHS1804/ARIHS1805/ ARIHS1807)
- · AnyOne (AOIHS2004/AOIHS2005/ AOIHS2007) · AnyRidge Octa 1(AROHS1604/
- AROHS1605/ AROHS1607)
- · Can get scan data without removing Scan Healing Abutment from Scan Post
- Different colors depend on the cuff size
- Scan healing abutment should be exposed 2.0mm on the surgical site for accurate scanning
- Scan Healing Abutment should be exposed 2.0mm from the surgical site for accurate scanning. Scanning would be much easier if you connect Scan Post when scanning seems difficult due to less exposure of Scan Healing Abutment or other conditions.
- · Select Scan Post based on the diameter of Scan Healing Abutment
- Scan Post is a disposable product and sold separately in batch of 10EA. for each package

## ARISH4004T SP4007.MTN 5	System		Profile Diameter	Scan Post	Height (mm)	Ref.C
## ARISH4007T ## ARISH5004T ## ARISH5004T ## ARISH5004T ## ARISH5004T ## ARISH5007T ## ARISH5007T ## ARISH6004T ## ARISH6004T ## ARISH6004T ## ARISH6005T ## ARISH6005T ## ARISH7004T ## ARISH7005T ## ARISH7005T ## ARISH7007T ## ARISH6004T ##						ARISH4004T
## ARISH5004T ## ARISH5004T ## ARISH5005T ## ARISH5005T ## ARISH5005T ## ARISH5007T ## ARISH6004T ## ARISH6004T ## ARISH6004T ## ARISH6004T ## ARISH6004T ## ARISH6005T ## ARISH6005T ## ARISH7004T ## ARISH7004T ## ARISH7004T ## ARISH7004T ## ARISH7007T ## ARISH7007T ## ARISH6004T ##			Ø4.0	SP4007.MTN	5	ARISH4005T
## AnyRidge ## Acceptable					7	ARISH4007T
AnyRidge Of .0 Cota 1 AnyRidge AnyRidge Of .0 Cota 1 AnyRidge AnyRidge Of .0 Cota 1 AnyRidge Of .0 AnyRidge Of .0 AnyRidge Of .0					4	ARISH5004T
AnyRidge AnyRidge AnyRidge AnyRidge AnyRidge AnyRidge AnyRidge Office of the part o			Ø5.0	SP5007.MTN	5	ARISH5005T
AnyRidge AnyRidge AnyRidge AnyRidge					7	ARISH5007T
AnyRidge					4	ARISH6004T
AnyRidge 07.0 SP7007.MTN 5 ARISH7004T 7 ARISH7005T 7 ARISH7007T 7 ARISH7007T 7 ARISH5004T 7 ARISH5004T 7 ARISH5005T 7 ARISH6005T 7 ARISH6004T 7 ARISH6004T 7 ARISH6005T 7 ARISH80005T 7 ARISH80005T 7 ARIOISHN4005T 7 ARIOISHN5005T 7 ARIOISHN5005T 7 ARIOISH8005T 7 ARIOISH8005T 7 ARIOISH8005T 7 ARIOISH8005T 7 ARIOISH85005T 7 ARIOISH85005T 7 ARIOISH85005T 7 ARIOISH85005T 7 ARIOISH86005T 7 ARIOISH86005T 7 ARIOISH86005T 7 ARIOISH86005T 7 ARIOISH86005T 7 ARIOISH86005T 7 ARIOISH87005T 7 ARIOISH87005			Ø6.0	SP6007.MTN	5	ARISH6005T
## ARISH70041 ARISH70041 ARISH7005T ARISH7007T ARISH7007T ARISH7007T ARISH7007T ARISH5004T ARISH5005T ARISH5005T ARISH5007T ARISH6004T ARISH6004T ARISH6004T ARISH6005T ARISH6005T ARISH6005T ARISH6007T ARISH6005T ARISH6007T ARISH6005T ARISHN4004T ARISHN5004T ARISHN5004T ARISHN5005T ARISHN5004T ARISHN5005T ARISHN5005T	. 5:				7	ARISH6007T
## ARISH7007T ## ARISH7007T ## ARISH7007T ## ARISH5004T ## ARNSH5005T ## ARNSH5005T ## ARNSH5007T ## ARNSH5007T ## ARNSH6004T ## ARNSH6004T ## ARNSH6004T ## AROISHN4004T ## AROISHN4004T ## AROISHN4004T ## AROISHN4004T ## AROISHN4004T ## AROISHR4004T ## AROISHR4004T ## AROISHR4004T ## AROISHR4004T ## AROISHR4004T ## AROISHR4004T ## AROISHR5004T ## AROISHR5004T ## AROISHR5004T ## AROISHR5004T ## AROISHR6004T ## AROISHR7004T ## A	AnyRic	ige			4	ARISH7004T
## ARNSH5004T SP5007.MTN 5			Ø7.0	SP7007.MTN	5	ARISH7005T
## ARNSH5005T Contract type SP5007.MTN 5					7	ARISH7007T
(Extra type) SP5007.MTN 5					4	ARNSH5004T
## ARNSH5007T ## ARNSH6004T ## ARNSH6004T ## ARNSH6004T ## ARNSH6005T ## ARNSH6005T ## ARNSH6005T ## ARNSH6005T ## ARNSH6005T ## ARNSH6007T ## ARNSH6007T ## AROISHN4004T ## AROISHN4005T ## AROISHN5004T ## AROISHN5004T ## AROISHN5005T ## AROISHN5005T ## AROISHN5005T ## AROISHN5005T ## AROISHR4004T ## AROISHR4005T ## AROISHR5004T ## AROISHR5005T ## AROISHR5005T ## AROISHR6005T ## AROISHR7004T ## AROISHR7005T ## AROISHR7				SP5007.MTN	5	ARNSH5005T
## AnyRidge Octa 1 Material Reservation			(Extra type)		7	ARNSH5007T
CExtra type SP6007.MTN 5					4	ARNSH6004T
AnyRidge				SP6007.MTN	5	ARNSH6005T
## AnyOne AnyOne AnyOne AnyOne AnyOne AnyOne AnyOne			(Extra type)			ARNSH6007T
## AROISHN4007T ## AROISHN4007T ## AROISHN5004T ## AROISHN5005T ## AROISHN5005T ## AROISHN5005T ## AROISHN5005T ## AROISHN5005T ## AROISHR4005T ## AROISHR4005T ## AROISHR4005T ## AROISHR4007T ## AROISHR5004T ## AROISHR5005T ## AROISHR5005T ## AROISHR5005T ## AROISHR6004T ## AROISHR6005T ## AROISHR6005T ## AROISHR6005T ## AROISHR6005T ## AROISHR7004T ## AROISHR7005T ## AROISHR4004T ## AOISH4005T ## AOISH4504T ## AOISH4505T ## AOISH4505T ## AROISHR5005T ## AROISHR5005T ## AROISHR7005T ## AROISHR7005T ## AROISHR7005T ## AROISHR7005T ## AROISHR4005T ## AROISHR4504T ## AOISH4505T ## AROISHR5005T					4	AROISHN4004T
## AROISHN5004T AROISHN5004T TAROISHN5004T TAROISHN5005T TAROISHN5005T TAROISHN5005T TAROISHN5005T TAROISHN5005T TAROISHN5005T TAROISHN5004T TAROISHN5004T TAROISHN5005T TAROISHN7004T TAROISHN5005T TAROISHN5005T		NC	Ø4.0	SP4007.MTN	5	AROISHN4005T
## AROISHN5004T AROISHN5004T AROISHN5005T AROISHN5005T AROISHN5007T AROISHN5007T AROISHN4004T AROISHN4005T AROISHN4005T AROISHN4007T AROISHN5004T AROISHN5004T AROISHN5005T AROISHN5005T AROISHN5005T AROISHN5005T AROISHN5005T AROISHN6004T AROISHN6005T AROISHN6005T AROISHN6005T AROISHN6005T AROISHN7004T AROISHN7005T AROISHN7005T AROISHN7005T AROISHN6005T ARO					7	AROISHN4007T
AnyRidge Octa 1 Rec			Ø5.0		4	AROISHN5004T
AnyRidge Octa 1 Rec				SP5007.MTN	5	AROISHN5005T
AnyRidge Octa 1 Representation of the province of the provinc					7	AROISHN5007T
AnyRidge Octa 1 Representation Octa 1 Oct					4	AROISHR4004T
Octa 1 Octa 1		•	Ø4.0	SP4007.MTN	5	AROISHR4005T
Octa 1 Ø5.0 SP5007.MTN 5 AROISHR5004T Ø6.0 SP6007.MTN 5 AROISHR6004T Ø7.0 SP7007.MTN 5 AROISHR6005T 7 AROISHR6005T 7 AROISHR6005T 7 AROISHR6007T 4 AROISHR7004T 4 AROISHR7004T 5 AROISHR7005T 7 AROISHR7005T 7 AROISHR0005T 7 AROISHR0005T 7 AROISHR7005T 7 AROISHR0005T 7 AROISHR0005T 4 AOISH4004T Ø4.0 SP4007.MTN 5 AOISH4005T 7 AOISH4504T Ø4.5 SP5007.MTN 5 AOISH4505T 7 AOISH4507T	AnyRidge				7	AROISHR4007T
Page	, ,				4	AROISHR5004T
AROISHR6004T AROISHR6004T AROISHR6004T AROISHR6005T 7 AROISHR6005T AROISHR7004T AROISHR7004T AROISHR7005T 7 AROISHR7005T 7 AROISHR7007T AROISH4004T AOISH4005T 7 AOISH4007T AOISH4504T AOISH4505T 7 AOISH4505T 7 AOISH4507T AOISH4507			Ø5.0	SP5007.MTN	5	AROISHR5005T
Ø6.0 SP6007.MTN 4 AROISHR6004T 5 AROISHR6005T 7 AROISHR6007T 4 AROISHR7004T 4 AROISHR7005T 7 AROISHR7007T 4 AOISH4004T 4 AOISH4005T 7 AOISH4007T 4 AOISH4504T 4 AOISH4504T 6 AOISH4505T 7 AOISH4507T			Ø5.0		7	AROISHR5007T
Ø7.0 SP7007.MTN 5 AROISHR7004T Ø4.0 SP4007.MTN 5 AROISHR7005T 7 AROISHR7007T 4 AOISH4004T Ø4.0 SP4007.MTN 5 AOISH4005T 7 AOISH4504T 4 AOISH4504T Ø4.5 SP5007.MTN 5 AOISH4505T AnvOne 7 AOISH4507T		W			4	AROISHR6004T
Ø7.0 SP7007.MTN 4 AROISHR7004T 5 AROISHR7005T 7 AROISHR7007T 4 AOISH4004T 4 AOISH4005T 7 AOISH4005T 7 AOISH4007T 4 AOISH4504T 4 AOISH4504T Ø4.5 SP5007.MTN 5 AOISH4505T AnvOne 7 AOISH4507T			Ø6.0	SP6007.MTN	5	AROISHR6005T
Ø7.0 SP7007.MTN 5 AROISHR7005T 7 AROISHR7007T 4 AOISH4004T 8 AOISH4005T 7 AOISH4007T 9 AOISH4504T 4 AOISH4504T 8 AOISH4505T 7 AOISH4507T					7	AROISHR6007T
Ø4.0 SP4007.MTN 5 AOISH4004T 7 AOISH4005T AOISH4005T 7 AOISH4007T 4 AOISH4504T Ø4.5 SP5007.MTN 5 AOISH4505T AnyOne 7 AOISH4507T					4	AROISHR7004T
Ø4.0 SP4007.MTN 4 AOISH4004T 5 AOISH4005T 7 AOISH4007T 4 AOISH4504T 4 AOISH4504T Ø4.5 SP5007.MTN 5 AOISH4505T AnyOne 7 AOISH4507T			Ø7.0	SP7007.MTN	5	AROISHR7005T
Ø4.0 SP4007.MTN 5 AOISH4005T 7 AOISH4007T 4 AOISH4504T Ø4.5 SP5007.MTN 5 AOISH4505T AnvOne 7 AOISH4507T					7	AROISHR7007T
7 AOISH4007T 4 AOISH4504T Ø4.5 SP5007.MTN 5 AOISH4505T 7 AOISH4507T					4	AOISH4004T
Ø4.5 SP5007.MTN 4 AOISH4504T AnvOne 4 AOISH4505T 7 AOISH4507T				SP4007.MTN	5	AOISH4005T
Ø4.5 SP5007.MTN 5 AOISH4505T AnvOne 7 AOISH4507T	AnyOne				7	AOISH4007T
AnyOne 7 AOISH4507T					4	AOISH4504T
AnvOne			Ø4.5	SP5007.MTN	5	AOISH4505T
AnyUne 4 AOISH5504T					7	AOISH4507T
					4	AOISH5504T
Ø5.5 SP6007.MTN 5 AOISH5505T			Ø5.5	SP6007.MTN	5	AOISH5505T
7 AOISH5507T					7	AOISH5507T
4 AOISH6504T					4	AOISH6504T
Ø6.5 SP7007.MTN 5 AOISH6505T			Ø6.5	SP7007.MTN	5	AOISH6505T
7 AOISH6507T					7	AOISH6507T



* If Scan Healing Abutment is exposed more than 2.5mm, it may unstablize a fixture and results in fixture failure.



Scan Post Carrier

	System	Length	Ref.C
С	ommom	19	SPC16



RP Analog Option

RP Analog

- · For Chairside/ Labside
- Included spare Abutment Screw
- Surpporting Dental CAD
- 3 Shape
- Exocad

System	Profile Diameter	Length (mm)	Туре	Ref.C
AnyRidge	Ø4.0	9	-	CANIAR4009
AnyRidge	Ø3.3	10	NC	AROLAN
Octa 1	Ø4.1	10	RC	AROLAR
AnyOne	Ø4.0	9	Only Ø3.5	CAOIAS3509
Internal	04.0	9	-	CAOIAR4009
	Ø3.5		Small	CEXEAS3509
AnyOne External	Ø4.1	9	Regular	CEXEAR4109
Extorria	Ø5.0		Wide	CEXEAW5009
AnyOne OneStage	Ø4.8	9	Cuff 1.8	OSRA18
MiNi	Ø3.0	9	-	CMIIAN3009
ST	Ø3.7	0	Small	OSRA3709
51	Ø4.3	9	Regular	OSRA4309
	Ø3.8		Small	OCTARA4
Octa Level	Ø4.8	9	Regular	OCTARA5
	Ø5.8		Wide	OCTARA6
MUA Level (N-Type)	Ø4.8	9	-	MUALA



ZrGEN Abutment Option

ZrGEN Abutment

- Abutment Screw included. . AnyRidge (AANMSF) . AnyOne Internal (AS20)

- . AnyOne Exeternal(SCS160/ RCS200) . AnyOne Stage (

- . MiNi (MIAZ1410) . ST(OSGSAS3110/OSGSAS3210) . Octa Level(IRCS200)
- · AnyRidge Octa 1(AROAS16B/ AROAS16)
- · Titanium Base
- 1Set(=Abutment 10ea)
 - included spare Abutment Screw
- MiNi ZrGEN has special ZrGEN Screw
- · Supporting DentalCAD
- 3 Shape
- Exocad
- Dental Wing
- Different groove number depend on the post size

-P=4.5 ▶ groove number : 2ea -P=5 ► groove number : 3ea -P=6 ► groove number : 4ea -P=8 ► groove number : 6ea



Standard

System	Diameter	Cuff Height	Post Height	Туре	Ref.C
			4.5		AANIPR4015.MTN
		0.6	6		AANIPR4016.MTN
			8		AANIPR4018.MTN
			4.5		AANIPR4025.MTN
		1.5	6		AANIPR4026.MTN
			8	Hex	AANIPR4028.MTN
			4.5	пех	AANIPR4035.MTN
		3.0	6		AANIPR4036.MTN
			8		AANIPR4038.MTN
			4.5		AANIPR4045.MTN
		4.0	6		AANIPR4046.MTN
	Ø4.0		8		AANIPR4048.MTN
	04.0		4.5		AANIPR4015N.MTN
		0.6	6		AANIPR4016N.MTN
			8		AANIPR4018N.MTN
			4.5		AANIPR4025N.MTN
		1.5	6		AANIPR4026N.MTN
			8	Non-Hex	AANIPR4028N.MTN
			4.5	Non-Hex	AANIPR4035N.MTN
		3.0	6		AANIPR4036N.MTN
			8		AANIPR4038N.MTN
			4.5		AANIPR4045N.MTN
		4.0	6		AANIPR4046N.MTN
AnyRidge			8		AANIPR4048N.MTN
7 ti iyi ilage		0.6	4.5	Hex	AANIPR4515.MTN
			6		AANIPR4516.MTN
			8		AANIPR4518.MTN
			4.5		AANIPR4525.MTN
		1.5	6		AANIPR4526.MTN
			8		AANIPR4528.MTN
			4.5		AANIPR4535.MTN
		3.0	6		AANIPR4536.MTN
			8		AANIPR4538.MTN
			4.5		AANIPR4545.MTN
		4.0	6	-	AANIPR4546.MTN
	Ø4.5		8		AANIPR4548.MTN
			4.5	-	AANIPR4515N.MTN
		0.6	6	-	AANIPR4516N.MTN
			8	-	AANIPR4518N.MTN
			4.5	-	AANIPR4525N.MTN
	-	1.5	6	-	AANIPR4526N.MTN
			8	Non-Hex	AANIPR4528N.MTN
			4.5	Non-nex	AANIPR4535N.MTN
		3.0	6		AANIPR4536N.MTN
			8		AANIPR4538N.MTN
			4.5		AANIPR4545N.MTN
		4.0	6		AANIPR4546N.MTN
			8		AANIPR4548N.MTN

Syst	tem	Diameter	Cuff Height	Post Height	Туре	Ref.C
			0.6			AROZGN4015.MTN
			1.5	4.5		AROZGN4025.MTN
			3.0	4.5		AROZGN4035.MTN
			4.0			AROZGN4045.MTN
			0.6			AROZGN4016.MTN
	NC	Ø4.0	1.5	6.0		AROZGN4026.MTN
		204.0	3.0	0.0		AROZGN4036.MTN
			4.0			AROZGN4046.MTN
			0.6			AROZGN4018.MTN
			1.5	8.0		AROZGN4028.MTN
	, 0		3.0	0.0		AROZGN4038.MTN
AnyRidge			4.0		_	AROZGN4048.MTN
Octa 1			0.6			AROZGR4515.MTN
			1.5	4.5		AROZGR4525.MTN
			3.0	7.0		AROZGR4535.MTN
		4.0			AROZGR4545.MTN	
	RC		0.6			AROZGR4516.MTN
		Ø4.5	1.5	6.0		AROZGR4526.MTN
		21.0	3.0	0.0		AROZGR4536.MTN
			4.0			AROZGR4546.MTN
			0.6			AROZGR4518.MTN
			1.5	8.0		AROZGR4528.MTN
			3.0	0.0		AROZGR4538.MTN
			4. 0			AROZGR4548.MTN
				5		AMUAPR5515N.MTN
			0.8	6		AMUAPR5516N.MTN
				8		AMUAPR5518N.MTN
				5		AMUAPR5525N.MTN
			1.7	6		AMUAPR5526N.MTN
MUA Level	Ø5.5		8	N-Type	AMUAPR5528N.MTN	
	Ø5.5		5	(Nobel)	AMUAPR5535N.MTN	
		3.0	6		AMUAPR5536N.MTN	
			8		AMUAPR5538N.MTN	
				5		AMUAPR5545N.MTN
			4.0	6		AMUAPR5546N.MTN
				8		AMUAPR5548N.MTN

Standard

Otaridald										
System	Diameter	Cuff Height	Post Height	Туре	Ref.C					
			4.5		AAOIPR4015.MTN					
		0.6	6		AAOIPR4016.MTN					
			8		AAOIPR4018.MTN					
			4.5		AAOIPR4025.MTN					
		1.5	6		AAOIPR4026.MTN					
			8	Hex	AAOIPR4028.MTN					
			4.5	I IOX	AAOIPR4035.MTN					
		3.0	6		AAOIPR4036.MTN					
			8		AAOIPR4038.MTN					
			4.5		AAOIPR4045.MTN					
		4.0	6		AAOIPR4046.MTN					
	Ø4.0		8		AAOIPR4048.MTN					
			4.5		AAOIPR4015N.MTN					
		0.6	6		AAOIPR4016N.MTN					
			8		AAOIPR4018N.MTN					
			4.5	-	AAOIPR4025N.MTN					
		1.5	6		AAOIPR4026N.MTN					
			8	Non-Hex	AAOIPR4028N.MTN					
			4.5		AAOIPR4035N.MTN					
		3.0	6		AAOIPR4036N.MTN					
			8		AAOIPR4038N.MTN					
		4.0	4.5		AAOIPR4045N.MTN					
AO			6		AAOIPR4046N.MTN					
AnyOne			8		AAOIPR4048N.MTN					
Internal		0.6	4.5		AAOIPR4515.MTN					
			6		AAOIPR4516.MTN					
			8 4.5		AAOIPR4518.MTN					
		1.5			AAOIPR4525.MTN					
		1.5	6		AAOIPR4526.MTN					
			8 4.5	Hex	AAOIPR4528.MTN					
		3.0			AAOIPR4535.MTN					
		3.0	6 8		AAOIPR4536.MTN AAOIPR4538.MTN					
			4.5		AAOIPR4538.IVITN					
		4.0	6		AAOIPR4545.MTN					
		4.0	8		AAOIPR4548.MTN					
	Ø4.5		4.5		AAOIPR4515N.MTN					
		0.6	6		AAOIPR4516N.MTN					
		0.0	8		AAOIPR4518N.MTN					
			4.5		AAOIPR4525N.MTN					
		1.5	6		AAOIPR4526N.MTN					
		1.5	8		AAOIPR4528N.MTN					
			4.5	Non-Hex	AAOIPR4528N.MTN					
		3.0	6		AAOIPR4536N.MTN					
		5.0	8		AAOIPR4538N.MTN					
			4.5		AAOIPR4538N.MTN					
		4.0	6		AAOIPR4546N.MTN					
		4.0	8							
			Ö		AAOIPR4548N.MTN					

Syste	em	Diameter	Cuff Height	Post Height	Туре	Ref.C
				4.5		AEXEPS4015.MTN
			0.6	6		AEXEPS4016.MTN
				8		AEXEPS4018.MTN
				4.5		AEXEPS4025.MTN
			1.5	6		AEXEPS4026.MTN
		Ø4.2		8		AEXEPS4028.MTN
		04.2		4.5		AEXEPS4035.MTN
			3.0	6		AEXEPS4036.MTN
				8		AEXEPS4038.MTN
				4.5		AEXEPS4045.MTN
			4.0	6		AEXEPS4046.MTN
	Small			8	_	AEXEPS4048.MTN
	Orrical			4.5		AEXEPS4515.MTN
			0.6	6		AEXEPS4516.MTN
				8		AEXEPS4518.MTN
				4.5		AEXEPS4525.MTN
			1.5	6		AEXEPS4526.MTN
		Ø4.5		8		AEXEPS4528.MTN
				4.5		AEXEPS4535.MTN
			3.0	6		AEXEPS4536.MTN
				8	-	AEXEPS4538.MTN
				4.5		AEXEPS4545.MTN
			4.0	6		AEXEPS4546.MTN
AnyOne				8	Hex	AEXEPS4548.MTN
External				4.5		AEXEPR4515.MTN
			0.6	6		AEXEPR4516.MTN
				8		AEXEPR4518.MTN
				4.5		AEXEPR4525.MTN
			1.5	6		AEXEPR4526.MTN
	Regular	Ø4.5		8		AEXEPR4528.MTN
			0.0	4.5		AEXEPR4535.MTN
			3.0	6		AEXEPR4536.MTN
				8 4.5		AEXEPR4538.MTN AEXEPR4545.MTN
			4.0	6		AEXEPR4546.MTN
			4.0	8		AEXEPR4548.MTN
				4.5		AEXEPW5515.MTN
			0.6	6		AEXEPW5516.MTN
			0.0	8		AEXEPW5518.MTN
				4.5		AEXEPW5525.MTN
			1.5	6		AEXEPW5526.MTN
			1.0	8		AEXEPW5528.MTN
	Wide	Ø5.5		4.5		AEXEPW5535.MTN
			3.0	6		AEXEPW5536.MTN
			0.0	8		AEXEPW5538.MTN
				4.5		AEXEPW5545.MTN
			4.0	6		AEXEPW5546.MTN
				8		AEXEPW5548.MTN
				4.5		AEXIPR5015.MTN
			0.6	6		AEXIPR5016.MTN
				8		AEXIPR5018.MTN
				4.5		AEXIPR5025.MTN
			1.5	6		AEXIPR5026.MTN
AnyOne	Cuff		-	8		AEXIPR5028.MTN
OneStage		Ø4.8		4.5	Octa	AEXIPR5035.MTN
3-			3.0	6		AEXIPR5036.MTN
			-	8		AEXIPR5038.MTN
				4.5		AEXIPR5045.MTN
			4.0	6		AEXIPR5046.MTN
			4.0	8		AEXIPR5048.MTN

Standard

Sys	tem	Diameter	Cuff Height	Post Height	Туре	Ref.C	Sys	stem	Diameter	Cuff Height		Туре	Ref.C
M	liNi	Ø3.0	0.6	2.5	Hex	MIPN3013.MTN				0.0	5		AOCEPS5015.MTN
	1			2.5	Non-Hex	MIPN3013N.MTN				0.8	6 8		AOCEPS5016.MTN
			0.6	4.5		OSGSPA3111.MTN OSGSPA3112.MTN					5		AOCEPS5018.MTN AOCEPS5025.MTN
			0.0	8	-	OSGSPA3113.MTN				1.7	6		AOCEPS5026.MTN
				4.5		OSGSPA3121.MTN					8	0-4-	AOCEPS5028.MTN
			1.5	6		OSGSPA3122.MTN					5	Octa	AOCEPS5035.MTN
				8	Hex	OSGSPA3123.MTN				3.0	6		AOCEPS5036.MTN
				4.5	110%	OSGSPA3131.MTN					8		AOCEPS5038.MTN
			3.0	6 8	_	OSGSPA3132.MTN				4.0	6		AOCEPS5045.MTN
				4.5		OSGSPA3133.MTN OSGSPA3141.MTN				4.0	8		AOCEPS5046.MTN AOCEPS5048.MTN
			4.0	6		OSGSPA3142.MTN		Small	Ø5.0		5		ANOEPS5015.MTN
	0	Ø4.0		8		OSGSPA3143.MTN				0.8	6		ANOEPS5016.MTN
	Small	<i>1</i> 04.0		4.5		OSGSPA3111N.MTN					8		ANOEPS5018.MTN
			0.6	6		OSGSPA3112N.MTN					5		ANOEPS5025.MTN
				8	_	OSGSPA3113N.MTN				1.7	6		ANOEPS5026.MTN
			1.5	4.5		OSGSPA3121N.MTN OSGSPA3122N.MTN					8 5	Non-Octa	ANOEPS5028.MTN ANOEPS5035.MTN
			1.5	8		OSGSPA3123N.MTN				3.0	6		ANOEPS5035.WTN
				4.5	Non-Hex	OSGSPA3131N.MTN				0.0	8		ANOEPS5038.MTN
			3.0	6		OSGSPA3132N.MTN					5		ANOEPS5045.MTN
				8		OSGSPA3133N.MTN				4.0	6		ANOEPS5046.MTN
				4.5		OSGSPA3141N.MTN					8		ANOEPS5048.MTN
			4.0	6		OSGSPA3142N.MTN				0.8	5		AOCEPR5515.MTN
				8 4.5		OSGSPA3143N.MTN OSGSPA3211.MTN				0.8	6 8		AOCEPR5516.MTN AOCEPR5518.MTN
			0.6	6	OSGSPA3212.MTN					5		AOCEPR5525.MTN	
			0.0	8		OSGSPA4018.MTN				1.7	6		AOCEPR5526.MTN
				4.5		OSGSPA4025.MTN					8	Ooto	AOCEPR5528.MTN
			1.5	6		OSGSPA4026.MTN					5	Octa	AOCEPR5535.MTN
				8	Hex	OSGSPA4028.MTN				3.0	6		AOCEPR5536.MTN
			4.5	110%	OSGSPA4035.MTN					8		AOCEPR5538.MTN	
			3.0	6		OSGSPA4036.MTN				4.0	5 6		AOCEPR5545.MTN
				8 4.5		OSGSPA4038.MTN OSGSPA4045.MTN	Octa			4.0	8		AOCEPR5546.MTN AOCEPR5548.MTN
			4.0	6	Non-Hex	OSGSPA4046.MTN	Level	Regular	Ø5.5		5		ANOEPR5515.MTN
ST		040		8		OSGSPA4048.MTN				0.8	6		ANOEPR5516.MTN
51		Ø4.0		4.5		OSGSPA3211N.MTN					8		ANOEPR5518.MTN
			0.6	6		OSGSPA3212N.MTN					5		ANOEPR5525.MTN
				8		OSGSPA4018N.MTN				1.7	6		ANOEPR5526.MTN
			1.5	4.5		OSGSPA4025N.MTN OSGSPA4026N.MTN					8 5	Non-Octa	ANOEPR5528.MTN ANOEPR5535.MTN
			1.5	8		OSGSPA4028N.MTN				3.0	6		ANOEPR5536.MTN
				4.5		OSGSPA4035N.MTN				0.0	8		ANOEPR5538.MTN
			3.0	6		OSGSPA4036N.MTN					5	6	ANOEPR5545.MTN
				8		OSGSPA4038N.MTN				4.0			ANOEPR5546.MTN
			4.0	4.5		OSGSPA4045N.MTN					8		ANOEPR5548.MTN
			4.0	6 8	_	OSGSPA4046N.MTN				0.8	5 6		AOCEPW6515.MTN
	Regular			4.5		OSGSPA4048N.MTN OSGSPA4515.MTN				0.0	8		AOCEPW6516.MTN AOCEPW6518.MTN
			0.6	6		OSGSPA4516.MTN					5		AOCEPW6525.MTN
				8		OSGSPA4518.MTN				1.7	6		AOCEPW6526.MTN
				4.5		OSGSPA3221.MTN					8	Octa	AOCEPW6528.MTN
			1.5	6	-	OSGSPA3222.MTN					5	Ocia	AOCEPW6535.MTN
				8	Hex	OSGSPA4528.MTN				3.0	6		AOCEPW6536.MTN
			3.0	4.5		OSGSPA4535.MTN OSGSPA4536.MTN					8 5		AOCEPW6538.MTN AOCEPW6545.MTN
			5.0	8	-	OSGSPA4538.MTN				4.0	6		AOCEPW6546.MTN
				4.5		OSGSPA4545.MTN			~~-		8		AOCEPW6548.MTN
			4.0	6		OSGSPA4546.MTN		Wide	Ø6.5		5		ANOEPW6515.MTN
	Ø4.5		8		OSGSPA4548.MTN				0.8	6		ANOEPW6516.MTN	
			4.5		OSGSPA4515N.MTN					8		ANOEPW6518.MTN	
		0.6	6		OSGSPA4516N.MTN				17	5		ANOEPW6525.MTN	
			8 4.5		OSGSPA4518N.MTN OSGSPA3221N.MTN				1.7	6 8		ANOEPW6526.MTN ANOEPW6528.MTN	
			1.5	6		OSGSPA3221N.MTN					5	Non-Octa	ANOEPW6535.MTN
			1.0	8	N	OSGSPA4528N.MTN				3.0	6		ANOEPW6536.MTN
				4.5 6	INON-HEX	OSGSPA4535N.MTN					8.0 6		ANOEPW6538.MTN
			3.0			OSGSPA4536N.MTN					5	5	ANOEPW6545.MTN
				8	8 4.5 6	OSGSPA4538N.MTN				4.0 6		ANOEPW6546.MTN	
			4.0			OSGSPA4545N.MTN					8	8	ANOEPW6548.MTN
			4.0			OSGSPA4546N.MTN							



Extra

System	Fixture Core	Diameter	Cuff Height	Post Height	Туре	Ref.C
				4.5		ARZXN4515.MTN
			0.6	6		ARZXN4516.MTN
				8		ARZXN4518.MTN
				4.5		ARZXN4525.MTN
			1.5	6		ARZXN4526 .MTN
				8	Hex	ARZXN4528 .MTN
				4.5	TIEX	ARZXN4535.MTN
			3.0	6		ARZXN4536.MTN
				8		ARZXN4538.MTN
				4.5		ARZXN4545.MTN
			4.0	6		ARZXN4546 .MTN
	Core 3.3	Ø4.5		8		ARZXN4548 .MTN
			0.6	4.5		ARZXN4515N.MTN
			0.0	6 8		ARZXN4516N.MTN ARZXN4518N.MTN
				4.5		ARZXN4525N.MTN
			1.5	6		ARZXN4526N.MTN
			1.5	8		ARZXN4528N.MTN
				4.5	Non -Hex	ARZXN4535N.MTN
			3.0	6		ARZXN4536N.MTN
			0.0	8		ARZXN4538N.MTN
				4.5		ARZXN4545N.MTN
			4.0	6		ARZXN4546N.MTN
				8		ARZXN4548N.MTN
				4.5		ARZXM503815.MTN
			0.6	6		ARZXM503816.MTN
				8		ARZXM503818.MTN
				4.5		ARZXM503825.MTN
			1.5	1.5 6		ARZXM503826.MTN
				8	Hex	ARZXM503828.MTN
				4.5	пех	ARZXM503835.MTN
			3.0	6		ARZXM503836.MTN
				8		ARZXM503838.MTN
				4.5		ARZXM503845.MTN
			4.0	6		ARZXM503846.MTN
AnyRidge		Ø5.0		8		ARZXM503848.MTN
			0.6	4.5		ARZXM503815N.MTN
			0.6	6		ARZXM503816N.MTN
				8 4.5		ARZXM503818N.MTN ARZXM503825N.MTN
			1.5	6		ARZXM503826N.MTN
			1.5	8	Non -Hex	ARZXM503828N.MTN
				4.5		ARZXM503835N.MTN
			3.0	6		ARZXM503836N.MTN
				8		ARZXM503838N.MTN
				4.5		ARZXM503845N.MTN
			4.0	6		ARZXM503846N.MTN
	Coros			8		ARZXM503848N.MTN
	Core3.8			4.5		ARZXM553815.MTN
			0.6	6		ARZXM553816.MTN
				8		ARZXM553818.MTN
				4.5		ARZXM553825.MTN
			1.5	6		ARZXM553826.MTN
				8	Hex	ARZXM553828.MTN
				4.5	110%	ARZXM553835.MTN
			3.0	6		ARZXM553836.MTN
				8		ARZXM553838.MTN
			4.0	4.5		ARZXM553845.MTN
			4.0	6		ARZXM553846.MTN
		Ø5.5		8		ARZXM553848.MTN
			0.6	4.5		ARZXM553815N.MTN
			0.0	6 8		ARZXM553816N.MTN ARZXM553818N.MTN
		,		4.5		ARZXM553825N.MTN
			1.5	6		ARZXM553826N.MTN
			1.0	8		ARZXM553828N.MTN
				4.5	Non -Hex	ARZXM553835N.MTN
			3.0	6		ARZXM553836N.MTN
			0.0	8		ARZXM553838N.MTN
				4.5		ARZXM553845N.MTN
			4.0	6		ARZXM553846N.MTN
				8		ARZXM553848N.MTN

Extra

System	Fixture Core	Diameter	Cuff Height	Post Height	Туре	Ref.C	System	Fixture Core	Diameter	Cuff Height	Post Height	Туре	Ref.C			
				4.5		ARZXM5015.MTN					4.5		ARZXL5515.MTN			
			0.6	6		ARZXM5016.MTN				0.6	6		ARZXL5516 .MTN			
				8		ARZXM5018.MTN					8		ARZXL5518 .MTN			
				4.5		ARZXM5025.MTN					4.5		ARZXL5525.MTN			
			1.5	6		ARZXM5026.MTN				1.5	6		ARZXL5526 .MTN			
				8	□ov	ARZXM5028.MTN					8	Hex	ARZXL5528 .MTN			
				4.5	Hex	ARZXM5035.MTN					4.5	TIEX	ARZXL5535.MTN			
			3.0	6		ARZXM5036.MTN				3.0	6		ARZXL5536 .MTI			
				8		ARZXM5038.MTN					8		ARZXL5538 .MTI			
				4.5		ARZXM5045.MTN			Ø5.5		4.5		ARZXL5545.MTI			
		Ø5.0	4.0	6		ARZXM5046.MTN				4.0	6		ARZXL5546 .MTI			
				8		ARZXM5048.MTN					8		ARZXL5548 .MTI			
	0.0ھ	25.0		4.5		ARZXM5015N.MTN					4.5		ARZXL5515N.MT			
						0.6	6		ARZXM5016N.MTN				0.6	6		ARZXL5516N.MT
				8		ARZXM5018N.MTN					8		ARZXL5518N.MT			
				4.5		ARZXM5025N.MTN					4.5		ARZXL5525N.MT			
			1.5	6		ARZXM5026N.MTN				1.5	6		ARZXL5526N.MT			
			8	Non Hov	ARZXM5028N.MTN					8	Non -Hex	ARZXL5528N.MT				
				4.5	Non -Hex	ARZXM5035N.MTN				3.0	4.5	NOIT-HEX	ARZXL5535N.MT			
			3.0	6		ARZXM5036N.MTN					6		ARZXL5536N.M			
				8		ARZXM5038N.MTN	AnyRidge				8		ARZXL5538N.M			
				4.5		ARZXM5045N.MTN					4.5		ARZXL5545N.MT			
			4.0	6		ARZXM5046N.MTN				4.0	6		ARZXL5546N.MT			
AnyRidge	Core4.0			8		ARZXM5048N.MTN		Core 4.8			8		ARZXL5548N.MT			
aryr llage	00104.0	64.0		4.5	Hex	ARZXM5515.MTN		0010 4.0		0.6	4.5	-	ARZXL6015.MTI			
			0.6	6		ARZXM5516.MTN					6		ARZXL6016.MTI			
				8		ARZXM5518.MTN					8		ARZXL6018.MTI			
				4.5		ARZXM5525.MTN					4.5		ARZXL6025.MTI			
			1.5	6		ARZXM5526.MTN				1.5	6	Hex	ARZXL6026.MTI			
				8		ARZXM5528.MTN					8		ARZXL6028.MTI			
				4.5	I ICA	ARZXM5535.MTN					4.5		ARZXL6035.MTI			
			3.0	6		ARZXM5536.MTN				3.0	6		ARZXL6036.MTI			
				8		ARZXM5538.MTN					8		ARZXL6038.MTI			
				4.5		ARZXM5545.MTN					4.5		ARZXL6045.MTI			
			4.0	6		ARZXM5546.MTN				4.0	6		ARZXL6046.MTI			
		Ø5.5		8		ARZXM5548.MTN			Ø6.0		8		ARZXL6048.MTI			
		20.0		4.5		ARZXM5515N.MTN			20.0		4.5		ARZXL6015N.MT			
			0.6	6		ARZXM5516N.MTN				0.6	6		ARZXL6016N.MT			
				8		ARZXM5518N.MTN					8		ARZXL6018N.MT			
				4.5		ARZXM5525N.MTN					4.5		ARZXL6025N.MT			
			1.5	6		ARZXM5526N.MTN				1.5	6	Non -Hex	ARZXL6026N.MT			
				8	Non -Hev	ARZXM5528N.MTN				3.0	8		ARZXL6028N.MT			
				4.5	Non -Hex	ARZXM5535N.MTN					4.5		ARZXL6035N.MT			
			3.0	6		ARZXM5536N.MTN					6		ARZXL6036N.MT			
					8		ARZXM5538N.MTN					8		ARZXL6038N.MT		
				4.5							4.5	-	ARZXL6045N.MT			
				6		ARZXM5546N.MTN				4.0	6		ARZXL6046N.MT			
				8		ARZXM5548N.MTN					8		ARZXL6048N.MT			

- ZrGEN Abutment

- Ti-base for Sirona Cerec users → CEREC
 In in Lab CAD Software, compatible with Xive Library







C-Type

Sys	stem	Diameter	Cuff Height	Post Height	Post Size	Ref.C
			0.5			ARCS3405.MTN
		Ø3.9	1			ARCS3410.MTN
			2		Cmall	ARCS3420.MTN
			0.5		Small	ARCS3805.MTN
Ar	nyRidge	Ø4.3	1	4.7		ARCS3810.MTN
			2			ARCS3820.MTN
			0.5			ARCL4505.MTN
		Ø5.5	1		Large	ARCL4510.MTN
			2			ARCL4520.MTN
			0.5			AROCSN3405.MTI
		Ø3.9	1.0			AROCSN3410.MTI
	NC		2.0		Cmall	AROCSN3420.MTN
			0.5		Small	AROCSN3805.MTN
		Ø4.3	1.0			AROCSN3810.MTI
			2.0			AROCSN3820.MTI
			0.5			AROCSR3405.MT
AnyRidge Octa 1		Ø3.9	1.0	4.5		AROCSR3410.MT
Ocia i			2.0			AROCSR3420.MT
			0.5		Small	AROCSR3805.MT
	RC .	Ø4.3	1.0			AROCSR3810.MT
			2.0			AROCSR3820.MT
			0.5			AROCLR4505.MT
		Ø5.5	1.0		Large	AROCLR4510.MT
			2.0			AROCLR4520.MT
			0.5			AOCS3405.MTN
		Ø3.9	1			AOCS3410.MTN
			2			AOCS3420.MTN
			0.5		Small	AOCS3805.MTN
А	nyOne	Ø4.3	1			AOCS3810.MTN
			2			AOCS3820.MTN
			0.5			AOCL4505.MTN
		Ø5.5	1		Large	AOCL4510.MTN
			2			AOCL4520.MTN
			0.5			STCSS3405.MTN
		Ø3.9	1			STCSS3410.MTN
	S		2			STCSS3420.MTN
	connection		0.5	4.7	Small	STCSS3805.MTN
		Ø4.3	1			STCSS3810.MTN
ST		2			STCSS3820.MTN	
			0.5			STCSR3405.MTN
		Ø3.9	1			STCSR3410.MTN
			2		_	STCSR3420.MTN
			0.5		Small	STCSR3805.MTN
	R	Ø4.3	1			STCSR3810.MTN
	connection		2			STCSR3820.MTN
			0.5			STCLR4505.MTN
		Ø5.5	1		Large	STCLR4510.MTN
		20.0			Large	C. JEI 11010.WITH

Scan Abutmet (C-type)

- Abutment Screw included. . AnyRidge (AANMSF) . AnyOne (AS20)


- . AnyRidge Octa 1(AROAS16B/ AROAS16)
- Scan Post for Sirona Cerec users CEREC
 In in Lab CAD Software, compatible with
 Xive Library

System		Profile Diameter	Cuff Height	Post Size	Ref.C
			0.5		ARICSS3405T
		Ø3.9	1		ARICSS3410T
			2		ARICSS3420T
			0.5	Small	ARICSS3805T
AnyRid	ge	Ø4.3	1		ARICSS3810T
			2		ARICSS3820T
			0.5		ARICSL4505T
		Ø5.5	1	Large	ARICSL4510T
			2		ARICSL4520T
			0.5		AROCSS3405NT
		Ø3.9	1		AROCSS3410NT
			2	0	AROCSS3420NT
	NC		0.5	Small	AROCSS3805NT
		Ø4.3	1		AROCSS3810NT
			2		AROCSS3820NT
	RC		0.5	Small	AROCSS3405RT
AnyRidge Octa 1		Ø3.9	1		AROCSS3410RT
Ooka i			2		AROCSS3420RT
			0.5		AROCSS3805RT
		Ø4.3	1		AROCSS3810RT
			2		AROCSS3820RT
			0.5		AROCSL4505RT
		Ø5.5	1	Large	AROCSL4510RT
			2		AROCSL4520RT
			0.5		AOICSS3405T
		Ø3.9	1		AOICSS3410T
			2	Small	AOICSS3420T
			0.5	SITIAII	AOICSS3805T
AnyOr	ne	Ø4.3	1		AOICSS3810T
			2		AOICSS3820T
			0.5		AOICSL4505T
		Ø5.5	1	Large	AOICSL4510T
			2		AOICSL4520T





○ TiGEN Prosthesis



○ TiGEN Abutment Option

TiGEN Abutment

- Abutment Screw included.
- . AnyRidge (AANMSF)
- . AnyOne Internal (AS20)
- . AnyOne Internal (AS20) . AnyOne Exeternal(SCS160/ RCS200) . AnyOne Stage (. MiNi (MIAZ1410) . ST(OSGSAS3110/OSGSAS3210)

- Octa Level(IRCS200)
- · AnyRidge Octa 1(AROAS16B/ AROAS16)
- · Pre-milled Abutment
- 1Set(=Abutment 10ea)
- included spare Abutment Screw
- Supporting DentalCAD
- 3Shape
- Exocad
- Dental Wings

Standard

System		Color Diameter Length		Туре	Ref.C	
AnyRidge		Gold	Ø10		Hex	ARTR1020.MTN
					Non-Hex	ARTR1020N.MTN
			Ø12		Hex	ARTR1220.MTN
					Non-Hex	ARTR1220N.MTN
AnyRidge Octa 1	NC	Gold	Ø10			AROTGN1020.MTN
			Ø12			AROTGN1220.MTN
	RC	Silver	Ø10			AROTGR1020.MTN
			Ø12			AROTGR1220.MTN
			Ø10		Hex	AOTR1020.MTN
Any	One	Pink	010		Non-Hex	AOTR1020N.MTN
Inte	nal		Ø12		Hex	AOTR1220.MTN
			012		Non-Hex	AOTR1220N.MTN
Anv/	AnyOne					AETS1220.MTN
Exte		N/A	Ø12		Hex	AETR1220.MTN
LXIC	irica					AETW1220.MTN
Mi	MiNi		Ø10		Hex	MITN1020.MTN
IVII					Non-Hex	MITN1020N.MTN
	Small		Ø10	20	Hex	OSTG3112.MTN
ST		Sky			Non-Hex	OSTG3112N.MTN
			Ø12		Hex	OSTG3111.MTN
					Non-Hex	OSTG3111N.MTN
01	Regular		Ø10		Hex	OSTG3212.MTN
					Non-Hex	OSTG3212N.MTN
			Ø12		Hex	OSTG3211.MTN
					Non-Hex	OSTG3211N.MTN
	Small	N/A	Ø10		Octa	OCTS1020.MTN
					Non-Octa	NOTS1020.MTN
			Ø12		Octa	OCTS1220.MTN
					Non-Octa	NOTS1220.MTN
	Regular		Ø10		Octa	OCTR1020.MTN
Octa					Non-Octa	NOTR1020.MTN
Level			Ø12		Octa	OCTR1220.MTN
					Non-Octa	NOTR1220.MTN
	Wide		Ø10		Octa	OCTW1020.MTN
					Non-Octa	NOTW1020.MTN
			Ø12		Octa	OCTW1220.MTN
					Non-Octa	NOTW1220.MTN



Extra EZ Connection

System	Color	Fixture Core	Diameter	Length	Туре	Ref.C
	Gold	3.3	Ø10	20	Hex	ARTXN1020.MTN
					Non-Hex	ARTXN1020N.MTN
			Ø12		Hex	ARTXN1220.MTN
					Non-Hex	ARTXN1220N.MTN
		4.0	Ø10		Hex	ARTXM1020.MTN
Any Didgo					Non-Hex	ARTXM1020N.MTN
AnyRidge			Ø12		Hex	ARTXM1220.MTN
					Non-Hex	ARTXM1220N.MTN
		4.8	Ø10		Hex	ARTXL1020.MTN
					Non-Hex	ARTXL1020N.MTN
			Ø12		Hex	ARTXL1220.MTN
					Non-Hex	ARTXL1220N.MTN





Why will AnyRidge work in any ridge?

Narran upper &

To maximize preservation by minimizing stress on the cortical bone.

**PEED®

- For faster, stronger osseointegration.
- New surface technology incorporating Ca²⁺ ions on the SLA treated surface.
- 100% elimination of any remaining acid from the conventional SLA process.

Wider fixture in a narkow crest

To maximize long term survival of implants.



- For smooth insertion and stronger primary stability.
- · No cuting edge for minimum invasion.
- · Ideal for soft bone cases.

Tapeled body

Excellent for simple installation

Excellent for simple installation and Immediate loading.

Narrow apical

For easier fixture insertion into a narrow ridge split incision

Case1









Case2





AnyRidge Clinical Case

Clinical Case 1

- Courtesy of Dr. Kwang-Bum Park

AnyRidge implant has excellent surface treatment which can be osseointegrated at this extreme case of bone defect.

Fig 1. This patient was 56 years old male patient and had a chief complaint of chewing difficulty on the left first mandibular molar due to periodontitis. On the panoramic radiograph, the tooth was floated with complete bone loss to the apex, and there was not enough bone to get initial stability for implant placement at the apex above mandibular nerve. So I decided to extract the tooth and wait for 4 months for regeneration of the socket.

Fig 2. The patient came back to my office after 4 months. Healing appeared good enough clinically, but the panoramic view still showed large socket defect. In many cases like this, we can expect some amount of bone filled in the socket which can allow minimal initial stability for implant placement.

Fig 3. When the flap was opened, I was very embarrassed that bone regeneration did not occur in the socket. None of remaining bone could be used for implant fixation.

Fig 4. A widest AnyRidge implant 8.0mm was placed on the mesial wall of extraction socket, but there was no initial stability. This trial was quite heroic treatment, but there was no other option except this because he spent many hours for this surgery.

Fig 5. The mixture of Allograft (Mega-Oss) and Synthetic bone (Bone Plus) was placed into the remaining socket defect and a collagen membrane was covered. Then primary closure was made with incision releasement on the periosteum.















Fig 6. On the panoramic view after surgery, we could find that none of the fixture was engaged with remaining bone, so it had more than 1.6mm gap from the tip to the depth of knife threads. I worried about the bone filling and success of the osseointegration on this fixture.

Fig 7. However, I was surprised with the hard cortical bone regeneration over the cover screw when I did the second stage surgery with the Biolaser.

Fig 8. On the intraoral radiograph taken several weeks after second surgery, we could see the fully regenerated bone into the bottom of thread depth.

Fig 9. The patient came back to our office to get one more implant on the maxillary upper molar after 2 years from the first implant surgery. The regenerated bone was matured and showed very stable crestal bone on the intraoral radiograph.













Clinical Case 2


- Courtesy of Dr. Kwang-Bum Park

Advantage of fuse abutment with AnyRidge implant for immediate loading in upper fully edentulous case

- Fig 1. An 80-year-old female patient presented with discomfort related to her upper teeth. About 10 years previously she underwent implant surgery in the mandible and received temporary teeth immediately as the bone density was sufficient for immediate loading. The patient requested a treatment plan for the upper arch that would give her immediate teeth.
- **Fig 2.** Clinical photos before surgery. The patient had no discomfort or complaints related to her mandibular implants. Plus her hygiene control was very good for maintaining healthy peri-implant tissue.
- **Fig 3.** All the remaining teeth were extracted. As shown, some teeth had severe periodontitis and some had decay at the cervix of the tooth. Drilling up to 2.9 mm was conducted at each implant site and eight 3.5×15mm implant fixtures were placed using a minimal flap design. All the fixtures showed excellent initial fixation, and the immediately-placed implants only had small socket defects.
- Fig 4. Eight fuse abutments were connected and the flaps were sutured to create a tight sealing against the fuse abutments
- **Fig 5.** The fuse abutments were prepared using a high speed handpiece for a temporary bridge that was already made before the implant surgery
- **Fig 6.** Panoramic scan taken immediately after surgery. The first premolar implant showed some mis-fit between the crown and the ratio.
- **Fig 7.** Intraoral scans taken 2 months after surgery. Shadows of the extraction socket can still be seen, but regenerated bone has started to fill the socket defects. The fuse abutments are functioning well without any problem. An impression was taken for customized zirconia abutments



















180203

Fig 12

Fig 8. Panoramic & intraoral scans taken 3 months after surgery. Zirconia customized abutments were connected to each fixture. The socket defects are completely filled with regenerated bone even in the case of immediate loading on the immediately placed implants.

Fig 9. Clinical photos of zirconia customized abutments and PMMA temporary bridge made using CAD/CAM technique. A zirconia abutment is excellent for both esthetics and hygiene maintenance. It has less than 1/10 bacterial accumulation on the surface compared with metals including titanium. PMMA provisional bridge is stronger than tooth resin, especially at the margin, so much beneficial for functional and occlusal tests.

Fig 10. A full zirconia one-piece bridge was made and delivered. The patient was very satisfied with the results, and thankful that she was provided with 'teeth' from the beginning to the end.

Fig 11. Panoramic scan of final restorations on day of delivery

Fig 12. Panoramic scans at 7 years follow-up

Clinical Case 3

- Courtesy of Dr. Soheil Bechara

Simultanious sinus lift and implant placement

Fig 1. The patient presented with huge bony defects around residual roots in the upper jaw. The treatment plan was to perform immediate implant placement and extract all decayed roots during one surgical session, as the patient had only one week to stay in the country.


Fig 2. Lateral window sinus lifting was performed on the right side with simultaneous implant placement in tooth area 17. In area 14 we can observe a huge bony defect which was thoroughly debrided until the margins of healthy bone. The Osteotomy was prepared with a 2mm final drill to place a 3.5x15mm Anyridge implant.


Fig 3. 3.5x15mm Anyridge implant was placed having only 2mm contact with the bone in the coronal part achieving 25 Ncm of initial torque.

Fig 4. The bony defect was filled with a Xenograft and covered with a collagen membrane. Although the defect is huge but it is still considered as an intra-bony defect with a good potential of bone regeneration.

Fig 5. On the left side two Anyridge implants were placed, immediate implant placement in area 24. Sinus lift with simultaneous implant placement in area 27.

Fig 6. No marginal bone loss, successful aesthetic and functional outcome.





Clinical Case 4

- Courtesy of Prof. Giuseppe Luongo

Immediate post-extraction insertion of implant and immediate loading.

Before Surgery





Fracture of #21 tooth. Good stability of the hard and soft tissues suggests immediate post-extraction insertion of implant and immediate loading.

Surgery









To protect the esthetic outcome of the procedure, the implant site was prepared via slightly palatal alveolar access.

A 4.5×11mm Anyridge was placed in the prepared site.






Biomaterial was added to the vestibular aspect to improve the stability of the esthetic outcome.

Temporary Prosthesis





12 Weeks after surgery



Tissue was ready to proceed with final abutment and crown.





The implant position was in harmony with the surrounding tissue and a prosthodontist completed the case using a zirconia framework.















Fig 2. Intraoral scan.

Fig 3, 4. Clinical photos of 1st surgery.

Fig 5, 6. Clinical photos of implant positioning.

Fig 7, 8. Clinical photo & scan of surgery.

Fig 9, 10. Clinical photos of immediate temporary crown in place.

Fig 11, 12. Clinical photos of healing and final abutment in place.

Fig 13. Clinical photo of zirconia framework in place.

Fig 14, 15. Clinical photo & intraoral scan of final crown at time of placement.

Fig 16, 17. Clinical photo & intraoral scan of final crown at 1-yr follow-up.

Fig 18, 19. Clinical photo & intraoral scan of final crown at 5-yr follow-up.





